

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Acura	MDX	2007 - 2013	7	Right	220 lbs	29.6	n/a	
Acura	MDX	2014 - 2020	7	Right	220 lbs	27.9	n/a	
Acura	RDX	2019 - Present	5	Right	440 lbs	27.5	24.8	
Audi	Q3	2015 - 2018	5	Left	175 lbs	n/a	n/a	MHFKR 264
Audi	Q3	2019 - Present	5	Right	220 lbs	24.3	n/a	MHFKR 427
Audi	Q5	2009 - 2017	5	Right	330 lbs	30.3	26.3	MHFKR 227
Audi	Q5	2018 - Present	5	Left	330 lbs	26.3	24.2	MHFKR 374
Audi	Q7	2016 - Present	7	Right	175 lbs	27.3	n/a	MHFKR 461
BMW	3 Series (Wagon)	2012 - 2018	5	Right	175 lbs	n/a	n/a	MHFKR 339
BMW	X1	2012 - 2015	5	Right	264 lbs	24	22	MHFKR 189
BMW	X1	2016 - Present	5	Right	175 lbs	25	n/a	MHFKR 342
BMW	X2	2019 - Present	5	Right	175 lbs	n/a	n/a	MHFK 439
BMW	X3	2004 - 2010	5	Right	440 lbs	30.5	28.5	MHFKR 208
BMW	X3	2011 - 2017	5	Right	440 lbs	29.2	27.1	
BMW	X3	2018 - Present	5	Left	175 lbs	25.3	n/a	MHFKR 448
BMW	X5	2019 - Present	5	Right	395 lbs	29.3	27.1	
Buick	Enclave	2008 - 2017	7	Right	440 lbs	28	25.9	
Buick	Enclave	2018 - Present	7	Right	440 lbs	29.1	27.3	
Buick	Encore	2013 - 2019	5	Right	264 lbs	26.8	24.8	MHFKR 252
Buick	Envision	2021 - Present	5	Left	175 lbs	23.8	n/a	
Buick	Rendezvous	2002 - 2007	5	Right	352 lbs	33.5	31.4	
Buick	Envision	2016 - 2020	5	Right	275 lbs	27.9	25.6	
Cadillac	SRX	2010 - 2016	5	Right	330 lbs	25.2	23.2	
Cadillac	XT5	2017 - Present	5	Right	220 lbs	27.1	n/a	
Cadillac	XT6	2020 - Present	7	Left	330 lbs	27.7	25.1	
Chevrolet	Blazer	2019 - Present	5	Right	220 lbs	27.5	n/a	
Chevrolet	Captiva	2007 - 2018	7	Right	220 lbs	30.4	n/a	
Chevrolet	City Express	2015 - 2018	7	Right	330 lbs	31.5	29.5	
Chevrolet	Equinox	2005 - 2009	5	Right	264 lbs	27.8	25.7	
Chevrolet	Equinox	2010 - 2017	5	Right	440 lbs	29.6	27.5	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit *universal option available
Chevrolet	Equinox	2018 - Present	5	Right	264 lbs	27.5	25	
Chevrolet	Express	1996 - Present	12	Right	352 lbs	33.5	31.4	
Chevrolet	HHR	2006 - 2011	5	Right	352 lbs	26.8	24.8	
Chevrolet	Orlando	2011 - 2018	7	Right	220 lbs	27.2	n/a	
Chevrolet	Suburban	2007 - 2014	7	Right	440 lbs	34	32	
Chevrolet	Suburban	2021 - Present	7	Right	275 lbs	29.5	28.3	
Chevrolet	Tahoe	2007 - 2014	7	Right	440 lbs	31.9	29.9	
Chevrolet	Tahoe	2015 - 2020	7	Right	440 lbs	27	25	
Chevrolet	Tahoe	2021 - Present	5	Right	220 lbs	29.9	n/a	
Chevrolet	Trailblazer	2002 - 2009	5	Right	264 lbs	28.8	26.7	
Chevrolet	Trailblazer	2021 - Present	5	Left	220 lbs	25.2	n/a	
Chevrolet	Traverse	2009 - 2017	7	Right	440 lbs	28	25.9	
Chevrolet	Traverse	2018 - Present	7	Right	440 lbs	31.4	29.1	
Chevrolet	Trax	2015 - Present	5	Right	264 lbs	26.8	24.8	
Chevrolet	Uplander	2005 - 2009	7	Right	440 lbs	35.1	33	
Chrysler	Pacifica	2017 - Present	7	Right	440 lbs	34.6	30.5	
Chrysler	Town & Country w/ out Stow n Go	2001 - 2007	7	Right	352 lbs	32.3	30.3	
Chrysler	Town & Country w/ Stow n Go	2001 - 2007	7	Right	440 lbs	30.3	28.3	
Chrysler	Town & Country w/ Stow n Go	2008 - Present	7	Right	440 lbs	35.6	33.8	
Dodge	Caravan & Grand Caravan w/ Stow n Go	2001 - 2007	7	Right	440 lbs	30.3	28.3	
Dodge	Caravan w/ out Stow n Go	2001 - 2007	7	Right	440 lbs	30.3	28.3	
Dodge	Durango	2004 - 2009	7	Right	330 lbs	29.6	27.5	
Dodge	Grand Caravan w/ Stow n Go	2008 - Present	7	Right	440 lbs	35.6	33.8	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Dodge	Journey	2009 - Present	5	Left	220 lbs	26.7	n/a	
Fiat	500L	2014 - Present	5	Right	330 lbs	28.4	26.3	
Fiat	500X	2016 - Present	5	Right	175 lbs	n/a	n/a	MHFKR 423
Ford	C-Max (gasoline version only)	2013 - 2018	5	Right	330 lbs	29.4	27.3	MHFKR 170
Ford	Ecosport	2018 - Present	5	Right	220 lbs	30.9	n/a	MHFKR 396
Ford	Edge	2015 - Present	5	Right	264 lbs	27.9	25.9	
Ford	Escape USA	2001 - 2007	5	Right	352 lbs	30.7	28.7	
Ford	Escape USA	2013 - 2019	5	Right	330 lbs	29.6	26.9	MHFKR 202
Ford	Escape USA	2020 - Present	5	Right	220 lbs	27.5	n/a	MHFKR 477
Ford	E-Series	1992 - 2014	12	Right	440 lbs	37.4	35.4	
Ford	Expedition	2007 - 2017	7	Right	352 lbs	30.3	28.3	
Ford	Expedition EL/Max	2007 - 2017	7	Right	440 lbs	31.1	29.1	
Ford	Explorer	2002 - 2005	5	Left	440 lbs	30.7	28.7	
Ford	Explorer	2011 - 2019	7	Right	330 lbs	27.9	25.9	
Ford	Explorer	2020 - Present	6	Left	440 lbs	30.5	27.5	
Ford	Flex	2013 - 2019	7	Right	440 lbs	29	27	
Ford	Focus Wagon	2012 - 2018	5	Right	330 lbs	28.2	26.1	
Ford	Freestar	2004 - 2007	7	Right	330 lbs	33.9	31.8	
Ford	Transit Connect LWB	2014 - 2019	7	Right	440 lbs	33.5	31.49	MHFKR 324
Ford	Transit Connect SWB	2014 - 2019	5	Right	440 lbs	33.5	31.4	
Genesis	GV80	2021 - Present	5	Right	220 lbs	26.7	n/a	
GMC	Acadia	2007 - 2017	7	Right	440 lbs	28	25.9	
GMC	Acadia	2017 - Present	7	Left	330 lbs	28.3	25.5	
GMC	Envoy	2002 - 2009	5	Right	264 lbs	28.8	26.7	
GMC	Savana	1996 - 2002	12	Right	352 lbs	33.5	31.4	
GMC	Terrain	2010 - 2017	5	Right	440 lbs	29.6	27.5	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
GMC	Terrain	2018 - Present	5	Right	264 lbs	27.5	25	
GMC	Yukon	2015 - 2020	7	Right	440 lbs	27	25	
GMC	Yukon	2021 - Present	5	Right	220 lbs	29.9	n/a	
GMC	Yukon Denali	2007 - 2014	7	Right	440 lbs	31.9	29.9	
GMC	Yukon XL	2007 - 2014	7	Right	440 lbs	34	32	
Honda	CR-V	2002 - 2006	5	Left	220 lbs		n/a	
Honda	CR-V	2007 - 2011	5	Right	220 lbs	31.5	n/a	
Honda	CR-V	2012 - 2016	5	Right	264 lbs	30.3	28.3	MHFKR 203
Honda	CR-V	2017 - Present	5	Right	440 lbs	30.9	28.7	
Honda	Fit	2015 - 2020	5	Right	330 lbs	27.6	25.5	MHFKR 321
Honda	HR-V	2016 - 2021	5	Right	440 lbs	26.6	24.6	MHFKR 313
Honda	Odyssey	2005 - 2010	7	Right	352 lbs	35.9	33.8	
Honda	Odyssey	2011 - 2013	7	Right	330 lbs	35.5	33.4	
Honda	Odyssey	2014 - 2017	7	Right	440 lbs	37	35	
Honda	Odyssey	2018 - Present	7	Right	440 lbs	33.4	31.4	
Honda	Passport	2019 - Present	5	Right	440 lbs	29.1	27.1	
Honda	Pilot	2009 - 2015	7	Right	440 lbs	29.6	27.5	
Honda	Pilot	2016 - Present	8	Right	330 lbs	29.1	26.3	
Honda	Fit	2009 - 2014	5	Right	220 lbs		n/a	
Hyundai	Entourage	2007 - 2009	7	Left	220 lbs		n/a	
Hyundai	Entourage w/ Stow n Go	2006 - 2009	7	Right	352 lbs	30.7	28.7	
Hyundai	Kona	2018 - Present	5	Right	175 lbs	25	n/a	MHFKR 431
Hyundai	Palisade	2020 - Present	7	Right	440 lbs	29.5	26.7	
Hyundai	Santa Fe	2001 - 2006	5	Right	220 lbs		n/a	
Hyundai	Santa Fe	2007 - 2012	5	Right	220 lbs	31.5	n/a	
Hyundai	Santa Fe	2013 - 2019	7	Right	330 lbs	27.2	25.1	MHFKR 271
Hyundai	Santa Fe	2021 - Present	5	Right	275 lbs	27.5	26.7	
Hyundai	Santa Fe	2019 - 2020	5	Right	275 lbs	27.5	25.1	
Hyundai	Santa Fe Sport	2013 - 2018	5	Right	264 lbs	27.8	25.7	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Hyundai	Tucson	2005 -2009	5	Left	330 lbs	28.4	26.3	
Hyundai	Tucson	2010 -2015	5	Right	330 lbs			
Hyundai	Tucson	2016 - 2021	5	Left	220 lbs	26.7	n/a	
Hyundai	Tucson HEV	2022 - Present	5	Left	220 lbs	27.5	n/a	
Jaguar	E-Pace	2018 - Present	5	Right	175 lbs	24.6	n/a	MHFKR 393
Jaguar	F Pace	2016 - Present	5	Right	175 lbs	26.1	n/a	MHFKR 361
Jeep	Cherokee	2014 - Present	5	Right	220 lbs	25.7	n/a	
Jeep	Compass	2018 - Present	5	Right	175 lbs	24.6	n/a	MHFKR 395
Jeep	Grand Cherokee	2005 - 2010	5	Right	220 lbs		n/a	
Jeep	Grand Cherokee	2011 - 2020	5	Right	264 lbs	28	25.9	MHFKR 254
Jeep	Grand Cherokee L	2021 - Present	6	Right	352 lbs	30.3	27.2	
Jeep	Liberty	2008 - 2012	5	Right	264 lbs	26.6	24.6	
Jeep	Renegade	2015 - 2020	5	Right	264 lbs	27	25	
Kia	Borrego	2009 - 2011	7	Right	440 lbs	30.9	28	
Kia	Carnival	2022 - Present	7	Right	440 lbs	35.2	31.8	
Kia	Niro BEV	2023 - Present	5	Left	175 lbs	n/a	n/a	MHFKR 558
Kia	Niro HEV	2017 - 2022	5	Left	175 lbs	22.4	n/a	
Kia	Rondo/Carens	2006 - 2012	5 or 7	Right	220 lbs			
Kia	Rondo/Carens	2013 - 2017	5	Right	220 lbs	27.6	n/a	MHFKR 233
Kia	Sedona	2002 - 2005	7	Left	220 lbs		n/a	
Kia	Sedona NON Stow n Go	2006 - 2014	7	Left	220 lbs		n/a	
Kia	Sedona with Stow n Go	2006 - 2014	7	Right	352 lbs	30.7	28.7	
Kia	Seltos	2021 - Present	5	Right	220 lbs	25.9	n/a	
Kia	Sorento	2003 - 2009	5	Right	220 lbs	30	n/a	
Kia	Sorento	2011 - 2015	7	Right	330 lbs	28.4	26.3	
Kia	Sorento	2016 - 2020	7	Left	330 lbs	27.8	25.7	MHFKR 319
Kia	Sorento	2021 - Present	7	Left	440 lbs	27.6	25.4	MHFKR 537
Kia	Soul	2010 - 2013	5	Right	220 lbs	26	n/a	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit *universal option available
Kia	Soul	2014 - 2019	5	Right	220 lbs	25.8	n/a	MHFKR 331
Kia	Sportage	2011 - 2016	5	Right	175 lbs		n/a	MHFKR 192
Kia	Sportage	2017 - Present	5	Left	175 lbs	25.9	n/a	MHFKR 332
Kia	Telluride	2020 - Present	8	Right	440 lbs	29.5	27.9	
Kia	Sedona	2015 - 2021	7	Right	440 lbs	33	27.5	
Kia	Soul	2020 - Present	5	Right	220 lbs	25.9	n/a	
Land Rover	Discovery 2	1999 - 2003	5	Left	220 lbs		n/a	
Land Rover	Discovery 5	2017 - Present	7	Right	220 lbs	30.7	n/a	MHFKR 363
Land Rover	Discovery Sport	2014 - 2019	7	Right	440 lbs	28.6	26.5	MHFKR 347
Land Rover	Freelander	2000 - 2005	5	Left	220 lbs		n/a	
Land Rover	Freelander	2006 - 2014	5	Right	330 lbs	28.2	26.1	
Land Rover	Range Rover Evoque	2011 - 2018	5	Right	220 lbs	24.8	n/a	MHFKR 190
Land Rover	Range Rover Evoque	2019 - Present	5	Right	175 lbs	24.2	n/a	
Land Rover	Range Rover Sport	2005 - 2013	5	Right	220 lbs	28.6	n/a	
Land Rover	Range Rover Sport	2014 - Present	5	Right	264 lbs	28	25.9	MHFKR 249
Land Rover	Range Rover Velar	2018 - Present	5	Right	175 lbs	23	n/a	MHFKR 487
Lexus	NX	2015 - Present	5	Right	175 lbs	24.2	n/a	MHFKR 419
Lexus	RX	2016 - Present	5	Left	175 lbs	24	n/a	MHFKR 334
Lexus	RX450h	2010 - 2015	5	Right	264 lbs	25.4	24.4	MHFKR 250
Lincoln	Corsair	2020 - Present	5	Right	220 lbs	25.5	n/a	
Lincoln	MKX	2016 - 2018	5	Right	275 lbs	25.5	24.4	
Mazda	5	2012 - 2015	7	Right	440 lbs	29.4	27.3	
Mazda	6	2003 - 2008	5	Right	220 lbs		n/a	
Mazda	CX-30	2020 - Present	5	Right	175 lbs	22.6	n/a	MHFKR 475
Mazda	CX-5	2013 - 2016	5	Right	330 lbs	30.3	28.3	MHFKR 228

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Mazda	CX-5	2017 - Present	5	Right	352 lbs	27.5	25.1	
Mazda	CX-9	2007 - 2015	7	Right	330 lbs	29.6	27.5	
Mazda	MX-30 EV	2022 - Present	5	Right	175 lbs	n/a	n/a	MHFKR 564
Mercedes	B Class	2006 - 2011	5	Right	220 lbs	27.6	n/a	
Mercedes	B Class	2012 - 2018	5	Right	264 lbs	27.3	25.1	MHFKR 195
Mercedes	B Class	2019 - Present	5	Right	175 lbs	26.7	n/a	MHFKR 436
Mercedes	C Class	2008 - 2014	5	Right	220 lbs		n/a	
Mercedes	C Class Wagon	2015 - Present	5	Right	220 lbs	23.7	n/a	MHFKR 318
Mercedes	E Class Estate	2017 - Present	5	Right	264 lbs	24	23	MHFKR 353
Mercedes	GL Class	2007 - 2012	7	Right	264 lbs	30.7	28.7	
Mercedes	GLA Class	2020 - Present	5	Right	220 lbs	26.3	n/a	MHFKR 499
Mercedes	GLA Class 45 AMG	2014 - 2019	5	Right	175 lbs	20	n/a	
Mercedes	GLC Class	2016 - Present	5	Left	175 lbs	27.4	n/a	MHFKR 424
Mercedes	GLE Class	2015 - 2018	5	Right	264 lbs	29.1	28.3	MHFKR 356
Mercedes	GLE Class	2019 - Present	5	Right	352 lbs	29.3	26.9	MHFKR 437
Mercedes	GLS	2016 -2019	7	Right	440 lbs	30.5	28.3	MHFKR 382
Mercedes	GLS	2020 - Present	5	Left	440 lbs	30.3	27.1	
Mercedes	M Class	2001 - 2005	5	Right	220 lbs		n/a	
Mercedes	M Class	2006 - 2011	5	Right	330 lbs	31.1	29.1	MHFKR 188
Mercedes	M Class	2012 - 2015	5	Right	352 lbs	30.3	28.3	MHFKR 201
MINI	Countryman	2017 - Mid 2020	5	Right	175 lbs	25.7	n/a	MHFKR 381
Mitsubishi	Eclipse Cross	2018 - Present	5	Right	175 lbs	26.1	n/a	MHFKR 388
Mitsubishi	Outlander	2014 - 2021	7	Right	330 lbs	28	25.9	MHFKR 209
Mitsubishi	Outlander PHEV	2014 - 2021	5	Right	330 lbs	28	25.9	
Nissan	Cube	2009 - 2014	5	Left	220 lbs	31.9	n/a	
Nissan	Murano	2003 - 2007	5	Left	440 lbs	28.5	26.7	MHFKR 401
Nissan	NV200	2010 - Present	7	Right	330 lbs	31.5	29.5	
Nissan	Pathfinder	2005 - 2012	7	Right	220 lbs	31.5	n/a	
Nissan	Pathfinder	2013 - 2021	7	Right	440 lbs	29.9	27.9	
Nissan	Pathfinder	2022 - Present	7	Right	330 lbs	29.1	26.3	

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Nissan	Quest	2011 - 2017	7	Right	440 lbs	34	32.5	
Nissan	Rogue	2014 - 2020	5	Right	330 lbs	28.9	27.5	MHFKR 274
Nissan	Rogue	2014 - 2019	7	Right	330 lbs	28.8	26.7	MHFKR 279
Nissan	Rouge	2021 - Present	5	Right	275 lbs	29.1	27.5	
Nissan	Rouge Sport/Qashqai	2017 - 2022	5	Right	220 lbs	23.6	n/a	MHFKR 262
Nissan	Xterra	2005 - 2012	5	Right	440 lbs	31.9	29.9	
Nissan	Versa Note	2014 - 2019	5	Right	220 lbs	25.6	n/a	MHFKR 240
Porche	Cayenne	2011 - 2018	5	Right	275 lbs	27.6	25.6	MHFKR 255
Ram	Promaster City	2015 - Present	5 & 7	Right	330 lbs	30	27.9	
Subaru	Ascent	2019 - Present	7	Right	275 lbs	29.3	26.7	
Subaru	Forester	2003 - 2008	5	Right	220 lbs	28	n/a	
Subaru	Forester	2014 - 2018	5	Left	330 lbs	29.6	27.5	MHFKR 333
Subaru	Forester	2019 - Present	5	Left	275 lbs	29.7	27.9	
Subaru	Forester XT	2012 - 2018	5	Left	330 lbs	29.2	27.1	
Subaru	Outback	2005 - 2009	5	Right	220 lbs	26.4	n/a	
Subaru	Outback	2010 - 2014	5	Left	220 lbs	27.6	n/a	MHFKR 327
Subaru	Outback	2015 - 2019	5	Right	330 lbs	27.5	25.5	
Subaru	Outback	2020 - Present	5	Left	220 lbs	27.5	n/a	
Tesla	Model X	2016 - Present	6	Right	175 lbs	23.8	n/a	
Tesla	Model Y	2019 - Present	5	Right	175 lbs	n/a	n/a	
Toyota	4 Runner	2010 - Present	5	Left	220 lbs	30.7		
Toyota	C-HR	2018 - Present	5	Left	175 lbs	n/a	n/a	MHFKR 379
Toyota	Highlander	2001 - 2007	5	Right	330 lbs	28.4	26.3	
Toyota	Highlander	2014 - 2019	7	Left	264 lbs	29.6	27.5	
Toyota	Highlander	2020 - Present	7	Left	220 lbs	27.9	n/a	
Toyota	Land Cruiser	1999 - 2007	8	Left	220 lbs		n/a	
Toyota	Land Cruiser	2008 - Present	7	Left	220 lbs	30.7	n/a	MHFKR 200
Toyota	Rav4	2013 - 2018	5	Right	352 lbs	28.4	26.3	MHFKR 241

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit <small>*universal option available</small>
Toyota	Rav4 (USA Version)	2006 - 2012	5	Left	440 lbs	34.3	32.2	MHFKR 211
Toyota	Rav4 (USA Version)	2019 - Present	5	Left	440 lbs	31.1	28.3	
Toyota	Sienna	2004 - 2010	7	Right	352 lbs	35.1	33	
Toyota	Sienna	2011 - 2020	7	Right	352 lbs	34.3	32.2	
Toyota	Sienna	2021 - Present	7	Right	440 lbs	34	31.5	
Toyota	Yaris (Hatchback)	2012 - 2018	5	Right	175 lbs	n/a	n/a	MHFKR 345
Volvo	V50	2004 - 2011	5	Right	220 lbs	23.7	n/a	
Volvo	V70	2001 - 2007	5	Right	220 lbs		n/a	
Volvo	V70	2008 - 2016	5	Right	264 lbs			
Volvo	V90	2017 - Present	5	Right	175 lbs	25.3	n/a	MHFKR 443
Volvo	XC40	2018 - Present	5	Left	175 lbs	26.2	n/a	MHFKR 414
Volvo	XC60	2009 - 2017	5	Right	352 lbs	27.2	25.1	
Volvo	XC60	2018 - Present	5	Left	175 lbs	28.3	n/a	MHFKR 386
Volvo	XC70	2001 - 2007	5	Right	220 lbs		n/a	
Volvo	XC70	2008 - 2016	5	Right	264 lbs			
Volvo	XC90	2015 - Present	7	Right	440 lbs	27.5	25.9	
VW	Atlas	2018 - Present	5	Right	275 lbs	30.3	27.7	
VW	Atlas Cross Sport	2020 - Present	5	Right	220 lbs	26.1	n/a	
VW	Golf Hatchback	2015 - 2020	5	Left	175 lbs	n/a	n/a	MHFKR 354
VW	Golf Mk8(eTSi)	2022 - Present	5	Left	175 lbs	n/a	n/a	MHFKR 559
VW	Golf Wagon	2015 - Present	5	Right	220 lbs	24.4	n/a	MHFKR 256
VW	Golf Wagon	2010 - 2014	5	Right	330 lbs	27.2	25.1	
VW	Taos	2022 - Present	5	Left	275 lbs	28	25.6	
VW	Tiguan (LWB)	2018 - Present	7	Right	220 lbs	25.9	n/a	MHFKR 422
VW	Touareg	2011 - 2018	5	Right	352 lbs	28.4	26.3	MHFKR 284