

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LM (inches)	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit *universal option available
Acura	MDX	2007 - 2013	7	Right	220 lbs		29.6	n/a	
Acura	MDX	2014 - 2020	7	Right	220 lbs		27.9	n/a	
Acura	MDX	2022 - Present	7	Right	220 lbs	29.5	27.5	n/a	
Acura	RDX	2019 - Present	5	Right	440 lbs		27.5	24.8	
Audi	Q3	2015 - 2018	5	Left	175 lbs	24.8	n/a	n/a	MHFKR 264
Audi	Q3	2019 - Present	5	Right	220 lbs	24.3	24.3	n/a	MHFKR 427
Audi	Q5	2009 - 2017	5	Right	330 lbs		30.3	26.3	MHFKR 227
Audi	Q5	2018 - Present	5	Left	330 lbs		26.3	24.2	MHFKR 374
Audi	Q7	2016 - Present	7	Right	175 lbs		27.3	n/a	MHFKR 461
BMW	3 Series (Wagon)	2012 - 2018	5	Right	175 lbs	23.8	n/a	n/a	MHFKR 339
BMW	X1	2012 - 2015	5	Right	264 lbs		24	22	MHFKR 189
BMW	X1	2016 - Present	5	Right	175 lbs		25	n/a	MHFKR 342
BMW	X2	2019 - Present	5	Right	175 lbs	26.5	n/a	n/a	MHFK 439
BMW	X3	2004 - 2010	5	Right	440 lbs		30.5	28.5	MHFKR 208
BMW	X3	2011 - 2017	5	Right	440 lbs		29.2	27.1	
BMW	X3	2018 - Present	5	Left	175 lbs	28.9	25.3	n/a	MHFKR 448
BMW	X5	2019 - Present	5	Right	395 lbs		29.3	27.1	
Buick	Enclave	2008 - 2017	7	Right	440 lbs		28	25.9	
Buick	Enclave	2018 - Present	7	Right	440 lbs		29.1	27.3	
Buick	Encore	2013 - 2019	5	Right	264 lbs		26.8	24.8	MHFKR 252
Buick	Encore GX	2020 - Present	5	Left	175 lbs	25.9	n/a	n/a	
Buick	Envision	2021 - Present	5	Left	175 lbs		23.8	n/a	
Buick	Rendezvous	2002 - 2007	5	Right	352 lbs		33.5	31.4	
Buick	Envision	2016 - 2020	5	Right	275 lbs		27.9	25.6	
Cadillac	SRX	2010 - 2016	5	Right	330 lbs		25.2	23.2	
Cadillac	XT5	2017 - Present	5	Right	220 lbs		27.1	n/a	
Cadillac	XT6	2020 - Present	7	Left	330 lbs		27.7	25.1	
Chevrolet	Blazer	2019 - Present	5	Right	220 lbs		27.5	n/a	
Chevrolet	Bolt EUV	2022 - Present	5	Left	175 lbs	n/a	22.8	n/a	
Chevrolet	Captiva	2007 - 2018	7	Right	220 lbs		30.4	n/a	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit *universal option available</b>
Chevrolet	City Express	2015 - 2018	7	Right	330 lbs		31.5	29.5	
Chevrolet	Equinox	2005 - 2009	5	Right	264 lbs		27.8	25.7	
Chevrolet	Equinox	2010 - 2017	5	Right	440 lbs		29.6	27.5	
Chevrolet	Equinox	2018 - Present	5	Right	264 lbs		27.5	25	
Chevrolet	Express	1996 - Present	12	Right	352 lbs		33.5	31.4	
Chevrolet	HHR	2006 - 2011	5	Right	352 lbs		26.8	24.8	
Chevrolet	Orlando	2011 - 2018	7	Right	220 lbs		27.2	n/a	
Chevrolet	Suburban	2007 - 2014	7	Right	440 lbs		34	32	
Chevrolet	Suburban	2021 - Present	7	Right	275 lbs		29.5	28.3	
Chevrolet	Tahoe	2007 - 2014	7	Right	440 lbs		31.9	29.9	
Chevrolet	Tahoe	2015 - 2020	7	Right	440 lbs		27	25	
Chevrolet	Tahoe	2021 - Present	5	Right	220 lbs		29.9	n/a	
Chevrolet	Trailblazer	2002 - 2009	5	Right	264 lbs		28.8	26.7	
Chevrolet	Trailblazer	2021 - Present	5	Left	220 lbs	26.7	25.2	n/a	
Chevrolet	Traverse	2009 - 2017	7	Right	440 lbs		28	25.9	
Chevrolet	Traverse	2018 - Present	7	Right	440 lbs		31.4	29.1	
Chevrolet	Trax	2015 - 2023	5	Right	264 lbs		26.8	24.8	
Chevrolet	Uplander	2005 - 2009	7	Right	440 lbs		35.1	33	
Chrysler	Pacifica	2017 - Present	7	Right	440 lbs		34.6	30.5	
Chrysler	Town & Country w/ out Stow n Go	2001 - 2007	7	Right	352 lbs		32.3	30.3	
Chrysler	Town & Country w/ Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	
Chrysler	Town & Country w/ Stow n Go	2008 - Present	7	Right	440 lbs		35.6	33.8	
Dodge	Caravan & Grand Caravan w/ Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	
Dodge	Caravan w/ out Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit *universal option available</b>
Dodge	Durango	2004 - 2009	7	Right	330 lbs		29.6	27.5	
Dodge	Grand Caravan w/ Stow n Go	2008 - Present	7	Right	440 lbs		35.6	33.8	
Dodge	Journey	2009 - Present	5	Left	220 lbs		26.7	n/a	
Dodge	Durango	2011 - Present	6	Right	275 lbs	29.5	29.5	26.9	
Fiat	500L	2014 - Present	5	Right	330 lbs		28.4	26.3	
Fiat	500X	2016 - Present	5	Right	175 lbs	21	n/a	n/a	MHFKR 423
Ford	Bronco Sport	2021 - Present	5	Right	220 lbs		26.5	n/a	
Ford	C-Max (gasoline version only)	2013 - 2018	5	Right	330 lbs		29.4	27.3	MHFKR 170
Ford	Ecosport	2018 - Present	5	Right	220 lbs		30.9	n/a	MHFKR 396
Ford	Edge	2007 - 2014	5	Right	440 lbs	31.5	28.3	26	
Ford	Edge	2015 - Present	5	Right	264 lbs		27.9	25.9	
Ford	Escape	2001 - 2007	5	Right	352 lbs		30.7	28.7	
Ford	Escape	2008 - 2012	5	Right	275 lbs		30.7	28.3	
Ford	Escape	2013 - 2019	5	Right	330 lbs		29.6	26.9	MHFKR 202
Ford	Escape	2020 - Present	5	Right	220 lbs		27.5	n/a	MHFKR 477
Ford	E-Series	1992 - 2014	12	Right	440 lbs		37.4	35.4	
Ford	Expedition	2007 - 2017	7	Right	352 lbs		30.3	28.3	
Ford	Expedition	2018 - Present	7	Right	220 lbs	34	29.9	n/a	
Ford	Expedition EL/Max	2007 - 2017	7	Right	440 lbs		31.1	29.1	
Ford	Explorer	2002 - 2005	5	Left	440 lbs		30.7	28.7	
Ford	Explorer	2011 - 2019	7	Right	330 lbs		27.9	25.9	
Ford	Explorer	2020 - Present	6	Left	440 lbs		30.5	27.5	
Ford	Flex	2013 - 2019	7	Right	440 lbs		29	27	
Ford	Focus Wagon	2012 - 2018	5	Right	330 lbs		28.2	26.1	
Ford	Freestar	2004 - 2007	7	Right	330 lbs		33.9	31.8	
Ford	Transit Connect LWB	2014 - 2023	7	Right	440 lbs		33.5	31.49	MHFKR 324

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit</b> *universal option available
Ford	Transit Connect SWB	2014 - 2023	5	Right	440 lbs		33.5	31.4	
Genesis	GV80	2021 - Present	5	Right	220 lbs		26.7	n/a	
GMC	Acadia	2007 - 2017	7	Right	440 lbs		28	25.9	
GMC	Acadia	2017 - Present	7	Left	330 lbs		28.3	25.5	
GMC	Envoy	2002 - 2009	5	Right	264 lbs		28.8	26.7	
GMC	Savana	1996 - 2002	12	Right	352 lbs		33.5	31.4	
GMC	Terrain	2010 - 2017	5	Right	440 lbs		29.6	27.5	
GMC	Terrain	2018 - Present	5	Right	264 lbs		27.5	25	
GMC	Yukon	2015 - 2020	7	Right	440 lbs		27	25	
GMC	Yukon	2021 - Present	5	Right	220 lbs		29.9	n/a	
GMC	Yukon Denali	2007 - 2014	7	Right	440 lbs		31.9	29.9	
GMC	Yukon XL	2007 - 2014	7	Right	440 lbs		34	32	
Honda	CR-V	2002 - 2006	5	Left	220 lbs			n/a	
Honda	CR-V	2007 - 2011	5	Right	220 lbs		31.5	n/a	
Honda	CR-V	2012 - 2016	5	Right	264 lbs		30.3	28.3	MHFKR 203
Honda	CR-V	2017 - 2022	5	Right	440 lbs		30.9	28.7	
Honda	CR-V	2023 - Present	5	Right	440 lbs	32.5	30.7	27.5	
Honda	Fit	2015 - 2020	5	Right	330 lbs		27.6	25.5	MHFKR 321
Honda	HR-V	2016 - 2022	5	Right	440 lbs		26.6	24.6	MHFKR 313
Honda	Odyssey	2005 - 2010	7	Right	352 lbs		35.9	33.8	
Honda	Odyssey	2011 - 2013	7	Right	330 lbs		35.5	33.4	
Honda	Odyssey	2014 - 2017	7	Right	440 lbs		37	35	
Honda	Odyssey	2018 - Present	7	Right	440 lbs		33.4	31.4	
Honda	Passport	2019 - Present	5	Right	440 lbs		29.1	27.1	
Honda	Pilot	2009 - 2015	7	Right	440 lbs		29.6	27.5	
Honda	Pilot	2016 - 2022	8	Right	330 lbs		29.1	26.3	
Honda	Pilot	2023 - Present	7	Right	275 lbs	29.1	28.3	25.9	
Honda	Fit	2009 - 2014	5	Right	220 lbs			n/a	
Hyundai	Entourage	2007 - 2009	7	Left	220 lbs			n/a	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit</b> *universal option available
Hyundai	Entourage w/ Stow n Go	2006 - 2009	7	Right	352 lbs		30.7	28.7	
Hyundai	Kona (Gas/HEV/BEV)	2018 - 2023	5	Right	175 lbs	27	25	n/a	MHFKR 431
Hyundai	Kona (mHEV)	2019 - 2023	5	Right	175 lbs	27	25	n/a	
Hyundai	Palisade	2020 - Present	7	Right	440 lbs		29.5	26.7	
Hyundai	Santa Fe	2001 - 2006	5	Right	220 lbs			n/a	
Hyundai	Santa Fe	2007 - 2012	5	Right	220 lbs		31.5	n/a	
Hyundai	Santa Fe	2019 - 2020	5	Right	275 lbs		27.5	25.1	
Hyundai	Santa Fe	2021 - Present	5	Right	275 lbs		27.5	26.7	
Hyundai	Santa Fe Sport	2013 - 2018	5	Right	264 lbs		27.8	25.7	
Hyundai	Tucson	2005 -2009	5	Left	330 lbs		28.4	26.3	
Hyundai	Tucson	2010 -2015	5	Right	330 lbs				
Hyundai	Tucson	2016 - 2021	5	Left	220 lbs		26.7	n/a	
Hyundai	Tucson HEV	2022 - Present	5	Left	220 lbs	29.5	27.5	n/a	
Jaguar	E-Pace	2018 - Present	5	Right	175 lbs		24.6	n/a	MHFKR 393
Jaguar	F Pace	2016 - Present	5	Right	175 lbs		26.1	n/a	MHFKR 361
Jeep	Cherokee	2014 - Present	5	Right	220 lbs		25.7	n/a	
Jeep	Compass	2018 - Present	5	Right	175 lbs		24.6	n/a	MHFKR 395
Jeep	Grand Cherokee	2005 - 2010	5	Right	220 lbs			n/a	
Jeep	Grand Cherokee	2011 - 2020	5	Right	264 lbs		28	25.9	MHFKR 254
Jeep	Grand Cherokee L	2021 - Present	6	Right	352 lbs		30.3	27.2	
Jeep	Liberty	2008 - 2012	5	Right	264 lbs		26.6	24.6	
Jeep	Patriot	2007 - 2017	5	Left	175 lbs	n/a	23.2	n/a	
Jeep	Renegade	2015 - 2020	5	Right	264 lbs		27	25	
Jeep	Grand Cherokee	2022 - Present	5	Right	330 lbs		28.9	26.5	
Jeep	Wagoneer/Grand Wagoneer	2022 - Present	8	Left	275 lbs		31.5	30.5	
Kia	Borrego	2009 - 2011	7	Right	440 lbs		30.9	28	
Kia	Carnival	2022 - Present	7	Right	440 lbs		35.2	31.8	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit *universal option available</b>
Kia	Niro BEV	2023 - Present	5	Left	175 lbs	25.9	n/a	n/a	MHFKR 558
Kia	Niro HEV	2017 - 2022	5	Left	175 lbs		22.4	n/a	
Kia	Rondo/Carens	2006 - 2012	5 or 7	Right	220 lbs				
Kia	Rondo/Carens	2013 - 2017	5	Right	220 lbs		27.6	n/a	MHFKR 233
Kia	Sedona	2002 - 2005	7	Left	220 lbs			n/a	
Kia	Sedona NON Stow n Go	2006 - 2014	7	Left	220 lbs			n/a	
Kia	Sedona with Stow n Go	2006 - 2014	7	Right	352 lbs		30.7	28.7	
Kia	Seltos	2021 - Present	5	Right	220 lbs		25.9	n/a	
Kia	Sorento	2003 - 2009	5	Right	220 lbs		30	n/a	
Kia	Sorento	2011 - 2015	7	Right	330 lbs		28.4	26.3	
Kia	Sorento	2016 - 2020	7	Left	330 lbs		27.8	25.7	MHFKR 319
Kia	Sorento	2021 - Present	7	Left	440 lbs		27.6	25.4	MHFKR 537
Kia	Soul	2010 - 2013	5	Right	220 lbs		26	n/a	
Kia	Soul	2014 - 2019	5	Right	220 lbs		25.8	n/a	MHFKR 331
Kia	Sportage	2011 - 2016	5	Right	175 lbs			n/a	MHFKR 192
Kia	Sportage	2017 - 2022	5	Left	175 lbs		25.9	n/a	MHFKR 332
Kia	Sportage LWB	2023 - Present	5	Left	275 lbs		29.1	27.1	
Kia	Telluride	2020 - Present	8	Right	440 lbs		29.5	27.9	
Kia	Sedona	2015 - 2021	7	Right	440 lbs		33	27.5	
Kia	Soul	2020 - Present	5	Right	220 lbs		25.9	n/a	
Land Rover	Discovery 2	1999 - 2003	5	Left	220 lbs			n/a	
Land Rover	Discovery 5	2017 - Present	7	Right	220 lbs		30.7	n/a	MHFKR 363
Land Rover	Discovery Sport	2014 - 2019	7	Right	440 lbs		28.6	26.5	MHFKR 347
Land Rover	Freelander	2000 - 2005	5	Left	220 lbs			n/a	
Land Rover	Freelander	2006 - 2014	5	Right	330 lbs		28.2	26.1	
Land Rover	Range Rover Evoque	2011 - 2018	5	Right	220 lbs		24.8	n/a	MHFKR 190

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit</b> *universal option available
Land Rover	Range Rover Evoque	2019 - Present	5	Right	175 lbs	29.1	24.2	n/a	
Land Rover	Range Rover Sport	2005 - 2013	5	Right	220 lbs		28.6	n/a	
Land Rover	Range Rover Sport	2014 - Present	5	Right	264 lbs		28	25.9	MHFKR 249
Land Rover	Range Rover Velar	2018 - Present	5	Right	175 lbs		23	n/a	MHFKR 487
Lexus	NX	2015 - 2021	5	Right	175 lbs		24.2	n/a	MHFKR 419
Lexus	RX	2010 - 2015	5	Right	264 lbs		25.4	24.4	MHFKR 250
Lexus	RX	2016 - 2022	5	Left	175 lbs		24	n/a	MHFKR 334
Lexus	RX	2023 - Present	5	Left	220 lbs	29.1	26.3	n/a	
Lexus	RX-L	2017 - 2022	5	Right	175 lbs	28	25.2	n/a	
Lincoln	Aviator	2020 - Present	7	Left	220 lbs		27.1	n/a	
Lincoln	Corsair	2020 - Present	5	Right	220 lbs		25.5	n/a	
Lincoln	MKT	2010 - 2019							
Lincoln	MKX	2016 - 2018	5	Right	275 lbs		25.5	24.4	
Mazda	5	2012 - 2015	7	Right	440 lbs		29.4	27.3	
Mazda	6	2003 - 2008	5	Right	220 lbs			n/a	
Mazda	CX-30	2020 - Present	5	Right	175 lbs		22.6	n/a	MHFKR 475
Mazda	CX-5	2013 - 2016	5	Right	330 lbs		30.3	28.3	MHFKR 228
Mazda	CX-5	2017 - Present	5	Right	352 lbs		27.5	25.1	
Mazda	CX-9	2007 - 2015	7	Right	330 lbs		29.6	27.5	
Mazda	MX-30 EV	2022 - Present	5	Right	175 lbs	25.6	n/a	n/a	MHFKR 564
Mercedes	B Class	2006 - 2011	5	Right	220 lbs		27.6	n/a	
Mercedes	B Class	2012 - 2018	5	Right	264 lbs		27.3	25.1	MHFKR 195
Mercedes	B Class	2019 - Present	5	Right	175 lbs		26.7	n/a	MHFKR 436
Mercedes	C Class	2008 - 2014	5	Right	220 lbs			n/a	
Mercedes	C Class Wagon	2015 - Present	5	Right	220 lbs		23.7	n/a	MHFKR 318
Mercedes	E Class Estate	2017 - Present	5	Right	264 lbs		24	23	MHFKR 353

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit *universal option available</b>
Mercedes	GL Class	2007 - 2012	7	Right	264 lbs		30.7	28.7	
Mercedes	GLA Class	2020 - Present	5	Right	220 lbs		26.3	n/a	MHFKR 499
Mercedes	GLA Class 45 AMG	2014 - 2019	5	Right	175 lbs		20	n/a	
Mercedes	GLB	2020 - Present	5	Right	220 lbs	30.1	27.9	n/a	MHFKR 584
Mercedes	GLC Class	2016 - Present	5	Left	175 lbs		27.4	n/a	MHFKR 424
Mercedes	GLE Class	2015 - 2018	5	Right	264 lbs		29.1	28.3	MHFKR 356
Mercedes	GLE Class	2019 - Present	5	Right	352 lbs		29.3	26.9	MHFKR 437
Mercedes	GLS	2016 -2019	7	Right	440 lbs		30.5	28.3	MHFKR 382
Mercedes	GLS	2020 - Present	5	Left	440 lbs		30.3	27.1	
Mercedes	M Class	2001 - 2005	5	Right	220 lbs			n/a	
Mercedes	M Class	2006 - 2011	5	Right	330 lbs		31.1	29.1	MHFKR 188
Mercedes	M Class	2012 - 2015	5	Right	352 lbs		30.3	28.3	MHFKR 201
Mercedes	Metris	2016 - 2023	8	Right	440 lbs		32.6	30.6	MHFKR 323
MINI	Countryman	2017 - Mid 2020	5	Right	175 lbs		25.7	n/a	MHFKR 381
Mitsubishi	Eclipse Cross	2018 - Present	5	Right	175 lbs		26.1	n/a	MHFKR 388
Mitsubishi	Outlander	2014 - 2021	7	Right	330 lbs		28	25.9	MHFKR 209
Mitsubishi	Outlander PHEV	2014 - 2021	5	Right	330 lbs		28	25.9	
Nissan	Cube	2009 - 2014	5	Left	220 lbs		31.9	n/a	
Nissan	Kicks	2018 - Present							
Nissan	Murano	2003 - 2007	5	Left	440 lbs		28.5	26.7	MHFKR 401
Nissan	NV200	2010 - Present	7	Right	330 lbs		31.5	29.5	
Nissan	Pathfinder	2005 - 2012	7	Right	220 lbs		31.5	n/a	
Nissan	Pathfinder	2013 - 2021	7	Right	440 lbs		29.9	27.9	
Nissan	Pathfinder	2022 - Present	7	Right	330 lbs		29.1	26.3	
Nissan	Quest	2011 - 2017	7	Right	440 lbs		34	32.5	
Nissan	Rogue	2014 - 2020	5	Right	330 lbs		28.9	27.5	MHFKR 274
Nissan	Rogue	2014 - 2019	7	Right	330 lbs		28.8	26.7	MHFKR 279
Nissan	Rouge	2021 - Present	5	Right	275 lbs		29.1	27.5	
Nissan	Rouge Sport/Qashqai	2017 - 2022	5	Right	220 lbs		23.6	n/a	MHFKR 262



**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

<b>Make</b>	<b>Model</b>	<b>Year</b>	<b>Seats</b>	<b>Fitting side</b>	<b>Max Device Weight</b>	<b>Max Device Height LM (inches)</b>	<b>Max Device Height LC (inches)</b>	<b>Max Device Height LP (inches)</b>	<b>Specific Restraint Kit</b> *universal option available
Nissan	Xterra	2005 - 2012	5	Right	440 lbs		31.9	29.9	
Nissan	Versa Note	2014 - 2019	5	Right	220 lbs		25.6	n/a	MHFKR 240
Porche	Cayenne	2011 - 2018	5	Right	275 lbs		27.6	25.6	MHFKR 255
Ram	Promaster City	2015 - Present	5 & 7	Right	330 lbs		30	27.9	
Subaru	Ascent	2019 - Present	7	Right	275 lbs		29.3	26.7	
Subaru	Crosstrek	2018 - 2023	5	Left	220 lbs	25.6	23.2		
Subaru	Forester	2003 - 2008	5	Right	220 lbs		28	n/a	
Subaru	Forester	2014 - 2018	5	Left	330 lbs		29.6	27.5	MHFKR 333
Subaru	Forester	2019 - Present	5	Left	275 lbs		29.7	27.9	
Subaru	Forester XT	2012 - 2018	5	Left	330 lbs		29.2	27.1	
Subaru	Outback	2005 - 2009	5	Right	220 lbs		26.4	n/a	
Subaru	Outback	2010 - 2014	5	Left	220 lbs		27.6	n/a	MHFKR 327
Subaru	Outback	2015 - 2019	5	Right	330 lbs		27.5	25.5	
Subaru	Outback	2020 - Present	5	Left	220 lbs		27.5	n/a	
Tesla	Model X	2016 - Present	6	Right	175 lbs		23.8	n/a	
Tesla	Model Y	2019 - Present	5	Right	175 lbs	25	n/a	n/a	
Toyota	4 Runner	2010 - Present	5	Left	220 lbs		30.7		
Toyota	C-HR	2018 - Present	5	Left	175 lbs	22.6	n/a	n/a	MHFKR 379
Toyota	Corolla Cross	2022 - Present	5	Left	175 lbs	27.2	24	n/a	
Toyota	Highlander	2001 - 2007	5	Right	330 lbs		28.4	26.3	
Toyota	Highlander	2014 - 2019	7	Left	264 lbs		29.6	27.5	
Toyota	Highlander	2020 - Present	7	Left	220 lbs		27.9	n/a	
Toyota	Grand Highlander	2024 - Present							
Toyota	Land Cruiser	1999 - 2007	8	Left	220 lbs			n/a	
Toyota	Land Cruiser	2008 - Present	7	Left	220 lbs		30.7	n/a	MHFKR 200
Toyota	Rav4	2013 - 2018	5	Right	352 lbs		28.4	26.3	MHFKR 241
Toyota	Rav4 (USA Version)	2006 - 2012	5	Left	440 lbs		34.3	32.2	MHFKR 211
Toyota	Rav4 (USA Version)	2019 - Present	5	Left	440 lbs		31.1	28.3	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits are available

Make	Model	Year	Seats	Fitting side	Max Device Weight	Max Device Height LM (inches)	Max Device Height LC (inches)	Max Device Height LP (inches)	Specific Restraint Kit *universal option available
Toyota	Sienna	2004 - 2010	7	Right	352 lbs		35.1	33	
Toyota	Sienna	2011 - 2020	7	Right	352 lbs		34.3	32.2	
Toyota	Sienna	2021 - Present	7	Right	440 lbs		34	31.5	
Toyota	Venza	2021 - Present							
Toyota	Yaris (Hatchback)	2012 - 2018	5	Right	175 lbs	25.6	n/a	n/a	MHFKR 345
Volvo	V50	2004 - 2011	5	Right	220 lbs		23.7	n/a	
Volvo	V70	2001 - 2007	5	Right	220 lbs			n/a	
Volvo	V70	2008 - 2016	5	Right	264 lbs				
Volvo	V90	2017 - Present	5	Right	175 lbs		25.3	n/a	MHFKR 443
Volvo	XC40	2018 - Present	5	Left	175 lbs		26.2	n/a	MHFKR 414
Volvo	XC60	2009 - 2017	5	Right	352 lbs		27.2	25.1	
Volvo	XC60	2018 - Present	5	Left	175 lbs		28.3	n/a	MHFKR 386
Volvo	XC70	2001 - 2007	5	Right	220 lbs			n/a	
Volvo	XC70	2008 - 2016	5	Right	264 lbs				
Volvo	XC90	2015 - Present	7	Right	440 lbs		27.5	25.9	
VW	Atlas	2018 - Present	5	Right	275 lbs		30.3	27.7	
VW	Atlas Cross Sport	2020 - Present	5	Right	220 lbs		26.1	n/a	
VW	Golf Hatchback	2015 - 2020	5	Left	175 lbs	23	n/a	n/a	MHFKR 354
VW	Golf Mk8(eTSi)	2022 - Present	5	Left	175 lbs	25.9	n/a	n/a	MHFKR 559
VW	Golf Wagon	2015 - Present	5	Right	220 lbs		24.4	n/a	MHFKR 256
VW	Golf Wagon	2010 - 2014	5	Right	330 lbs		27.2	25.1	
VW	ID.4	2021 - Present	5	Left	175 lbs	28.7	24.6		
VW	Taos	2022 - Present	5	Left	275 lbs		28	25.6	
VW	Tiguan (LWB)	2018 - Present	7	Right	220 lbs		25.9	n/a	MHFKR 422
VW	Tiguan (SWB)	2018 - Present	5	Left	220 lbs		27.3	n/a	MHFKR 359
VW	Touareg	2011 - 2018	5	Right	352 lbs		28.4	26.3	MHFKR 284