

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Acura | MDX | 2007 - 2013 | 7 | Right | 220 lbs | | 29.6 | N/A | |
| Acura | MDX | 2014 - 2020 | 7 | Right | 220 lbs | | 27.9 | N/A | |
| Acura | MDX | 2022 - 2026 | 7 | Right | 220 lbs | 29.5 | 27.5 | N/A | |
| Acura | RDX | 2019 - 2026 | 5 | Right | 440 lbs | | 27.5 | 24.8 | |
| Audi | Q3 | 2015 - 2018 | 5 | Left | 175 lbs | 24.8 | n/a | N/A | MHFKR 264 |
| Audi | Q3 | 2019 - 2025 | 5 | Right | 220 lbs | 24.3 | 24.3 | N/A | MHFKR 427 |
| Audi | Q4 E-Tron | 2022 - 2026 | 5 | Right | 175 lbs | 25.7 | 23 | N/A | |
| Audi | Q5 | 2009 - 2017 | 5 | Right | 330 lbs | | 30.3 | 26.3 | MHFKR 227 |
| Audi | Q5 | 2018 - 2025 | 5 | Left | 330 lbs | | 26.3 | 24.2 | MHFKR 374 |
| Audi | Q7 | 2016 - 2025 | 7 | Right | 175 lbs | | 27.3 | N/A | MHFKR 461 |
| BMW | 3 Series (Wagon) | 2012 - 2018 | 5 | Right | 175 lbs | 23.8 | n/a | N/A | MHFKR 339 |
| BMW | X1 | 2012 - 2015 | 5 | Right | 264 lbs | | 24 | 22 | MHFKR 189 |
| BMW | X1 | 2016 - 2022 | 5 | Right | 175 lbs | | 25 | N/A | MHFKR 342 |
| BMW | X2 | 2019 - 2025 | 5 | Right | 175 lbs | 26.5 | n/a | N/A | MHFK 439 |
| BMW | X3 | 2004 - 2010 | 5 | Right | 440 lbs | | 30.5 | 28.5 | MHFKR 208 |
| BMW | X3 | 2011 - 2017 | 5 | Right | 440 lbs | | 29.2 | 27.1 | |
| BMW | X3 | 2018 - 2024 | 5 | Left | 175 lbs | 28.9 | 25.3 | N/A | MHFKR 448 |
| BMW | X5 | 2019 - 2026 | 5 | Right | 395 lbs | | 29.3 | 27.1 | |
| BMW | X7 | 2019 - 2026 | 7 | Right | 220 lbs | N/A | 29.2 | N/A | |
| Buick | Enclave | 2008 - 2017 | 7 | Right | 440 lbs | | 28 | 25.9 | |
| Buick | Enclave | 2018 - 2024 | 7 | Right | 440 lbs | | 29.1 | 27.3 | |
| Buick | Enclave | 2025 - 2028 | 6 | Right | 220 lbs | N/A | 26.3 | N/A | |
| Buick | Encore | 2013 - 2022 | 5 | Right | 264 lbs | | 26.8 | 24.8 | MHFKR 252 |
| Buick | Encore GX | 2020 - 2026 | 5 | Left | 175 lbs | 25.9 | n/a | N/A | |
| Buick | Envision | 2021 - 2026 | 5 | Left | 175 lbs | | 23.8 | N/A | |
| Buick | Envista | 2024 - 2028 | 5 | Left | 175 lbs | N/A | 22.4 | N/A | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Buick | Rendezvous | 2002 - 2007 | 5 | Right | 352 lbs | | 33.5 | 31.4 | |
| Buick | Envision | 2016 - 2020 | 5 | Right | 275 lbs | | 27.9 | 25.6 | |
| Cadillac | Lyric | 2023 - 2027 | 5 | Right | 175 lbs | 28.8 | 25.3 | N/A | |
| Cadillac | SRX | 2010 - 2016 | 5 | Right | 330 lbs | | 25.2 | 23.2 | |
| Cadillac | XT5 | 2017 - 2026 | 5 | Right | 220 lbs | | 27.1 | N/A | |
| Cadillac | XT6 | 2020 - 2025 | 7 | Left | 330 lbs | | 27.7 | 25.1 | |
| Chevrolet | Blazer | 2019 - 2025 | 5 | Right | 220 lbs | | 27.5 | N/A | |
| Chevrolet | Blazer EV | 2024 - 2027 | 5 | Right | 175 lbs | N/A | 24 | N/A | |
| Chevrolet | Bolt EUV | 2022 - 2023 | 5 | Left | 175 lbs | n/a | 22.8 | N/A | |
| Chevrolet | Captiva | 2007 - 2018 | 7 | Right | 220 lbs | | 30.4 | N/A | |
| Chevrolet | City Express | 2015 - 2018 | 7 | Right | 330 lbs | | 31.5 | 29.5 | |
| Chevrolet | Equinox | 2005 - 2009 | 5 | Right | 264 lbs | | 27.8 | 25.7 | |
| Chevrolet | Equinox | 2010 - 2017 | 5 | Right | 440 lbs | | 29.6 | 27.5 | |
| Chevrolet | Equinox | 2018 - 2024 | 5 | Right | 264 lbs | | 27.5 | 25 | |
| Chevrolet | Equinox | 2025 - 2028 | 5 | Left | 220 lbs | N/A | 26.3 | N/A | |
| Chevrolet | Express | 1996 - 2026 | 12 | Right | 352 lbs | | 33.5 | 31.4 | |
| Chevrolet | HHR | 2006 - 2011 | 5 | Right | 352 lbs | | 26.8 | 24.8 | |
| Chevrolet | Orlando | 2011 - 2018 | 7 | Right | 220 lbs | | 27.2 | N/A | |
| Chevrolet | Suburban | 2007 - 2014 | 7 | Right | 440 lbs | | 34 | 32 | |
| Chevrolet | Suburban | 2021 - 2026 | 7 | Right | 275 lbs | | 29.5 | 28.3 | |
| Chevrolet | Tahoe | 2007 - 2014 | 7 | Right | 440 lbs | | 31.9 | 29.9 | |
| Chevrolet | Tahoe | 2015 - 2020 | 7 | Right | 440 lbs | | 27 | 25 | |
| Chevrolet | Tahoe | 2021 - 2025 | 5 | Right | 220 lbs | | 29.9 | N/A | |
| Chevrolet | Trailblazer | 2002 - 2009 | 5 | Right | 264 lbs | | 28.8 | 26.7 | |
| Chevrolet | Trailblazer | 2021 - 2026 | 5 | Left | 220 lbs | 26.7 | 25.2 | N/A | |
| Chevrolet | Traverse | 2009 - 2017 | 7 | Right | 440 lbs | | 28 | 25.9 | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|--------------------------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Chevrolet | Traverse | 2018 - 2023 | 7 | Right | 440 lbs | | 31.4 | 29.1 | |
| Chevrolet | Traverse | 2024 - Split Gen | 6-Jul | Right | 440 or 220 | N/A | 31.4 or 26.3 | 29.1 or N/A | |
| Chevrolet | Traverse | 2025 - 2028 | 6 | Right | 220 lbs | N/A | 26.3 | N/A | |
| Chevrolet | Trax | 2015 - 2023 | 5 | Right | 264 lbs | | 26.8 | 24.8 | |
| Chevrolet | Trax | 2024 - 2028 | 5 | Left | 175 lbs | N/A | 21.6 | N/A | |
| Chevrolet | Uplander | 2005 - 2009 | 7 | Right | 440 lbs | | 35.1 | 33 | |
| Chrysler | Pacifica | 2017 - 2026 | 7 | Right | 440 lbs | | 34.6 | 30.5 | |
| Chrysler | Town & Country w/ out Stow n Go | 2001 - 2007 | 7 | Right | 352 lbs | | 32.3 | 30.3 | |
| Chrysler | Town & Country w/ Stow n Go | 2001 - 2007 | 7 | Right | 440 lbs | | 30.3 | 28.3 | |
| Chrysler | Town & Country w/ Stow n Go | 2008 - 2016 | 7 | Right | 440 lbs | | 35.6 | 33.8 | |
| Dodge | Caravan & Grand Caravan w/ Stow n Go | 2001 - 2007 | 7 | Right | 440 lbs | | 30.3 | 28.3 | |
| Dodge | Caravan w/ out Stow n Go | 2001 - 2007 | 7 | Right | 440 lbs | | 30.3 | 28.3 | |
| Dodge | Durango | 2004 - 2009 | 7 | Right | 330 lbs | | 29.6 | 27.5 | |
| Dodge | Durango | 2011 - 2026 | 6 | Right | 275 lbs | 29.5 | 29.5 | 26.9 | |
| Dodge | Grand Caravan w/ Stow n Go | 2008 - 2020 | 7 | Right | 440 lbs | | 35.6 | 33.8 | |
| Dodge | Journey | 2009 -2020 | 5 | Left | 220 lbs | | 26.7 | N/A | |
| Fiat | 500L | 2014 - 2022 | 5 | Right | 330 lbs | | 28.4 | 26.3 | |
| Fiat | 500X | 2016 - 2022 | 5 | Right | 175 lbs | 21 | n/a | n/a | MHFKR 423 |
| Ford | Bronco Sport | 2021 - 2026 | 5 | Right | 220 lbs | | 26.5 | N/A | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|-------------------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Ford | C-Max (gasoline version only) | 2013 - 2018 | 5 | Right | 330 lbs | | 29.4 | 27.3 | MHFKR 170 |
| Ford | Ecosport | 2018 - 2022 | 5 | Right | 220 lbs | | 30.9 | N/A | MHFKR 396 |
| Ford | Edge | 2007 - 2014 | 5 | Right | 440 lbs | 31.5 | 28.3 | 26 | |
| Ford | Edge | 2015 - 2024 | 5 | Right | 264 lbs | | 27.9 | 25.9 | |
| Ford | Escape | 2001 - 2007 | 5 | Right | 352 lbs | | 30.7 | 28.7 | |
| Ford | Escape | 2008 - 2012 | 5 | Right | 275 lbs | | 30.7 | 28.3 | |
| Ford | Escape | 2013 - 2019 | 5 | Right | 330 lbs | | 29.6 | 26.9 | MHFKR 202 |
| Ford | Escape | 2020 - 2026 | 5 | Right | 220 lbs | | 27.5 | N/A | MHFKR 477 |
| Ford | E-Series | 1992 - 2014 | 12 | Right | 440 lbs | | 37.4 | 35.4 | |
| Ford | Expedition | 2007 - 2017 | 7 | Right | 352 lbs | | 30.3 | 28.3 | |
| Ford | Expedition | 2018 - 2024 | 7 | Right | 220 lbs | 34 | 29.9 | N/A | |
| Ford | Expedition EL/Max | 2007 - 2017 | 7 | Right | 440 lbs | | 31.1 | 29.1 | |
| Ford | Explorer | 2002 - 2005 | 5 | Left | 440 lbs | | 30.7 | 28.7 | |
| Ford | Explorer | 2011 - 2019 | 7 | Right | 330 lbs | | 27.9 | 25.9 | |
| Ford | Explorer | 2020 - 2026 | 6 | Left | 440 lbs | | 30.5 | 27.5 | |
| Ford | Flex | 2013 - 2019 | 7 | Right | 440 lbs | | 29 | 27 | |
| Ford | Focus Wagon | 2012 - 2018 | 5 | Right | 330 lbs | | 28.2 | 26.1 | |
| Ford | Freestar | 2004 - 2007 | 7 | Right | 330 lbs | | 33.9 | 31.8 | |
| Ford | Mustang Mach-E | 2021 - 2025 | 5 | Right | 175 lbs | 27.1 | 24 | | |
| Ford | Transit Connect LWB | 2014 - 2023 | 7 | Right | 440 lbs | | 33.5 | 31.49 | MHFKR 324 |
| Ford | Transit Connect SWB | 2014 - 2023 | 5 | Right | 440 lbs | | 33.5 | 31.4 | |
| Genesis | GV80 | 2021 - 2026 | 5 | Right | 220 lbs | | 26.7 | N/A | |
| GMC | Acadia | 2007 - 2017 | 7 | Right | 440 lbs | | 28 | 25.9 | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|--------------|---------------|--------------------|--------------------|---------------------------|----------------------------|---|--|--|---|
| GMC | Acadia | 2017 - 2023 | 7 | Left | 330 lbs | | 28.3 | 25.5 | |
| GMC | Acadia | 2024 - 2028 | 6 | Right | 220 lbs | N/A | 26.3 | N/A | |
| GMC | Envoy | 2002 - 2009 | 5 | Right | 264 lbs | | 28.8 | 26.7 | |
| GMC | Savana | 1996 - 2002 | 12 | Right | 352 lbs | | 33.5 | 31.4 | |
| GMC | Terrain | 2010 - 2017 | 5 | Right | 440 lbs | | 29.6 | 27.5 | |
| GMC | Terrain | 2018 - 2024 | 5 | Right | 264 lbs | | 27.5 | 25 | |
| GMC | Terrain | 2025 - 2027 | 5 | Left | 220 lbs | N/A | 26.3 | N/A | |
| GMC | Yukon | 2015 - 2020 | 7 | Right | 440 lbs | | 27 | 25 | |
| GMC | Yukon | 2021 - 2026 | 5 | Right | 220 lbs | | 29.9 | N/A | |
| GMC | Yukon Denali | 2007 - 2014 | 7 | Right | 440 lbs | | 31.9 | 29.9 | |
| GMC | Yukon XL | 2007 - 2014 | 7 | Right | 440 lbs | | 34 | 32 | |
| Honda | CR-V | 2002 - 2006 | 5 | Left | 220 lbs | | | N/A | |
| Honda | CR-V | 2007 - 2011 | 5 | Right | 220 lbs | | 31.5 | N/A | |
| Honda | CR-V | 2012 - 2016 | 5 | Right | 264 lbs | | 30.3 | 28.3 | MHFKR 203 |
| Honda | CR-V | 2017 - 2022 | 5 | Right | 440 lbs | | 30.9 | 28.7 | |
| Honda | CR-V | 2023 - 2026 | 5 | Right | 440 lbs | 31.1 | 29.8 | 27.6 | |
| Honda | Fit | 2015 - 2020 | 5 | Right | 220 lbs | | 27.6 | N/A | MHFKR 321 |
| Honda | HR-V | 2016 - 2022 | 5 | Right | 440 lbs | | 26.6 | 24.6 | MHFKR 313 |
| Honda | HR-V | 2023 - 2027 | 5 | Right | 220 lbs | 30.7 | 26.9 | | |
| Honda | Odyssey | 2005 - 2010 | 7 | Right | 352 lbs | | 35.9 | 33.8 | |
| Honda | Odyssey | 2011 - 2013 | 7 | Right | 330 lbs | | 35.5 | 33.4 | |
| Honda | Odyssey | 2014 - 2017 | 7 | Right | 440 lbs | | 37 | 35 | |
| Honda | Odyssey | 2018 - 2026 | 7 | Right | 440 lbs | | 33.4 | 31.4 | |
| Honda | Passport | 2019 - 2026 | 5 | Right | 440 lbs | | 29.1 | 27.1 | |
| Honda | Pilot | 2009 - 2015 | 7 | Right | 440 lbs | | 29.6 | 27.5 | |
| Honda | Pilot | 2016 - 2022 | 8 | Right | 330 lbs | | 29.1 | 26.3 | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|------------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Honda | Pilot | 2023 - 2026 | 7 | Right | 275 lbs | 29.1 | 28.3 | 25.9 | |
| Honda | Fit | 2009 - 2014 | 5 | Right | 220 lbs | | | N/A | |
| Hyundai | Entourage | 2007 - 2009 | 7 | Left | 220 lbs | | | N/A | |
| Hyundai | Entourage w/ Stow n Go | 2006 - 2009 | 7 | Right | 352 lbs | | 30.7 | 28.7 | |
| Hyundai | Ioniq 5 | 2022 - 2026 | 5 | Left | 175 lbs | 25.8 | N/A | N/A | |
| Hyundai | Kona (Gas/HEV/BEV) | 2018 - 2023 | 5 | Right | 175 lbs | 27 | 25 | N/A | MHFKR 431 |
| Hyundai | Kona (mHEV) | 2019 - 2023 | 5 | Right | 175 lbs | 27 | 25 | N/A | |
| Hyundai | Kona (Gas) | 2024 - 2026 | 5 | Right | 220 lbs | 28.3 | 26.5 | N/A | |
| Hyundai | Palisade | 2020 - 2025 | 7 | Right | 440 lbs | | 29.5 | 26.7 | |
| Hyundai | Santa Fe | 2001 - 2006 | 5 | Right | 220 lbs | | | N/A | |
| Hyundai | Santa Fe | 2007 - 2012 | 5 | Right | 220 lbs | | 31.5 | N/A | |
| Hyundai | Santa Fe | 2019 - 2020 | 5 | Right | 275 lbs | | 27.5 | 25.1 | |
| Hyundai | Santa Fe | 2021 - 2023 | 5 | Right | 275 lbs | | 27.5 | 26.7 | |
| Hyundai | Santa Fe | 2024 - 2028 | 7 | Left | 440 lbs | N/A | 29.3 | 27.3 | |
| Hyundai | Santa Fe Sport | 2013 - 2018 | 5 | Right | 264 lbs | | 27.8 | 25.7 | |
| Hyundai | Tucson | 2005 - 2009 | 5 | Left | 330 lbs | | 28.4 | 26.3 | |
| Hyundai | Tucson | 2010 - 2015 | 5 | Right | 330 lbs | | | | |
| Hyundai | Tucson | 2016 - 2021 | 5 | Left | 220 lbs | | 26.7 | N/A | |
| Hyundai | Tucson HEV | 2022 - 2026 | 5 | Left | 220 lbs | 29.5 | 27.5 | N/A | |
| Infiniti | QX60 | 2022 - 2026 | 8 | Right | 220 lbs | 30.0 | 28.0 | N/A | |
| Jaguar | E-Pace | 2018 - 2026 | 5 | Right | 175 lbs | | 24.6 | N/A | MHFKR 393 |
| Jaguar | F Pace | 2016 - 2026 | 5 | Right | 175 lbs | | 26.1 | N/A | MHFKR 361 |
| Jeep | Cherokee | 2014 - 2023 | 5 | Right | 220 lbs | | 25.7 | N/A | |
| Jeep | Compass | 2018 - 2025 | 5 | Right | 175 lbs | | 24.6 | N/A | MHFKR 395 |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|-------------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Jeep | Grand Cherokee | 2005 - 2010 | 5 | Right | 220 lbs | | | N/A | |
| Jeep | Grand Cherokee | 2011 - 2020 | 5 | Right | 264 lbs | | 28 | 25.9 | MHFKR 254 |
| Jeep | Grand Cherokee | 2022 - 2026 | 5 | Right | 330 lbs | | 28.9 | 26.5 | |
| Jeep | Grand Cherokee L | 2021 - 2026 | 6 | Right | 352 lbs | | 30.3 | 27.2 | |
| Jeep | Liberty | 2008 - 2012 | 5 | Right | 264 lbs | | 26.6 | 24.6 | |
| Jeep | Patriot | 2007 - 2017 | 5 | Left | 175 lbs | n/a | 23.2 | N/A | |
| Jeep | Renegade | 2015 - 2020 | 5 | Right | 264 lbs | | 27 | 25 | |
| Jeep | Wagoneer/Grand Wagoneer | 2022 - 2026 | 8 | Left | 275 lbs | | 31.5 | 30.5 | |
| Kia | Borrego | 2009 - 2011 | 7 | Right | 440 lbs | | 30.9 | 28 | |
| Kia | Carnival | 2022 - 2026 | 7 | Right | 440 lbs | | 35.2 | 31.8 | |
| Kia | Niro HEV | 2017 - 2022 | 5 | Left | 175 lbs | | 22.4 | N/A | |
| Kia | Rondo/Carens | 2006 - 2012 | 5 or 7 | Right | 220 lbs | | | | |
| Kia | Rondo/Carens | 2013 - 2017 | 5 | Right | 220 lbs | | 27.6 | N/A | MHFKR 233 |
| Kia | Sedona | 2002 - 2005 | 7 | Left | 220 lbs | | | N/A | |
| Kia | Sedona NON Stow n Go | 2006 - 2014 | 7 | Left | 220 lbs | | | N/A | |
| Kia | Sedona with Stow n Go | 2006 - 2014 | 7 | Right | 352 lbs | | 30.7 | 28.7 | |
| Kia | Seltos | 2021 - 2026 | 5 | Right | 220 lbs | | 25.9 | N/A | |
| Kia | Sorento | 2003 - 2009 | 5 | Right | 220 lbs | | 30 | N/A | |
| Kia | Sorento | 2011 - 2015 | 7 | Right | 330 lbs | | 28.4 | 26.3 | |
| Kia | Sorento | 2016 - 2020 | 7 | Left | 330 lbs | | 27.8 | 25.7 | MHFKR 319 |
| Kia | Sorento | 2021 - 2026 | 7 | Left | 440 lbs | | 27.6 | 25.4 | MHFKR 537 |
| Kia | Soul | 2010 - 2013 | 5 | Right | 220 lbs | | 26 | N/A | |
| Kia | Soul | 2014 - 2019 | 5 | Right | 220 lbs | | 25.8 | N/A | MHFKR 331 |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Kia | Soul | 2020 - 2026 | 5 | Right | 220 lbs | | 25.9 | N/A | |
| Kia | Sportage | 2011 - 2016 | 5 | Right | 175 lbs | | | N/A | MHFKR 192 |
| Kia | Sportage | 2017 - 2022 | 5 | Left | 175 lbs | | 25.9 | N/A | MHFKR 332 |
| Kia | Sportage LWB | 2023 - 2026 | 5 | Left | 275 lbs | | 29.1 | 27.1 | |
| Kia | Telluride | 2020 - 2026 | 8 | Right | 440 lbs | | 29.5 | 27.9 | |
| Kia | Sedona | 2015 - 2021 | 7 | Right | 440 lbs | | 33 | 27.5 | |
| Land Rover | Discovery 2 | 1999 - 2003 | 5 | Left | 220 lbs | | | N/A | |
| Land Rover | Discovery 5 | 2017 - 2026 | 7 | Right | 220 lbs | | 30.7 | N/A | MHFKR 363 |
| Land Rover | Discovery Sport | 2014 - 2019 | 7 | Right | 440 lbs | | 28.6 | 26.5 | MHFKR 347 |
| Land Rover | Freelander | 2000 - 2005 | 5 | Left | 220 lbs | | | N/A | |
| Land Rover | Freelander | 2006 - 2014 | 5 | Right | 330 lbs | | 28.2 | 26.1 | |
| Land Rover | Range Rover Evoque | 2011 - 2018 | 5 | Right | 220 lbs | | 24.8 | N/A | MHFKR 190 |
| Land Rover | Range Rover Evoque | 2019 - 2023 | 5 | Right | 175 lbs | 29.1 | 24.2 | N/A | |
| Land Rover | Range Rover Sport | 2005 - 2013 | 5 | Right | 220 lbs | | 28.6 | N/A | |
| Land Rover | Range Rover Sport | 2014 - 2022 | 5 | Right | 264 lbs | | 28 | 25.9 | MHFKR 249 |
| Land Rover | Range Rover Velar | 2018 - 2026 | 5 | Right | 175 lbs | | 23 | N/A | MHFKR 487 |
| Lexus | NX | 2015 - 2021 | 5 | Right | 175 lbs | | 24.2 | N/A | MHFKR 419 |
| Lexus | NX | 2022 - 2026 | 5 | Right | 220 lbs | 27.6 | 24.5 | N/A | |
| Lexus | RX | 2010 - 2015 | 5 | Right | 264 lbs | | 25.4 | 24.4 | MHFKR 250 |
| Lexus | RX | 2016 - 2022 | 5 | Left | 175 lbs | | 24 | N/A | MHFKR 334 |
| Lexus | RX | 2023 - 2026 | 5 | Left | 220 lbs | 29.1 | 26.3 | N/A | |
| Lexus | RX-L | 2017 - 2022 | 5 | Right | 175 lbs | 28 | 25.2 | N/A | |
| Lincoln | Aviator | 2020 - 2026 | 7 | Left | 220 lbs | | 27.1 | N/A | |
| Lincoln | Corsair | 2020 - 2026 | 5 | Right | 220 lbs | | 25.5 | N/A | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Lincoln | MKT | 2010 - 2019 | 6 | Right | 220 lbs | N/A | 27.5 | N/A | |
| Lincoln | MKX | 2016 - 2018 | 5 | Right | 275 lbs | | 25.5 | 24.4 | |
| Lincoln | Nautilus | 2019 - 2023 | 5 | Right | 220 lbs | 28.7 | 25.5 | N/A | |
| Lincoln | Nautilus | 2024 - 2026 | 5 | Right | 220 lbs | N/A | 25.9 | N/A | |
| Mazda | 5 | 2012 - 2015 | 7 | Right | 440 lbs | | 29.4 | 27.3 | |
| Mazda | 6 | 2003 - 2008 | 5 | Right | 220 lbs | | | N/A | |
| Mazda | CX-30 | 2020 - 2026 | 5 | Right | 175 lbs | | 22.6 | N/A | MHFKR 475 |
| Mazda | CX-50 | 2023 - 2026 | 5 | Right | 175 lbs | 28.7 | 26.3 | N/A | |
| Mazda | CX-70 | 2025 - 2028 | 5 | Right | 220 lbs | 29.5 | 26.7 | N/A | |
| Mazda | CX-5 | 2013 - 2016 | 5 | Right | 330 lbs | | 30.3 | 28.3 | MHFKR 228 |
| Mazda | CX-5 | 2017 - 2025 | 5 | Right | 352 lbs | | 27.5 | 25.1 | |
| Mazda | CX-9 | 2007 - 2015 | 7 | Right | 330 lbs | | 29.6 | 27.5 | |
| Mazda | MX-30 EV | 2022 - 2026 | 5 | Right | 175 lbs | 25.6 | n/a | N/A | MHFKR 564 |
| Mercedes | B Class | 2006 - 2011 | 5 | Right | 220 lbs | | 27.6 | N/A | |
| Mercedes | B Class | 2012 - 2018 | 5 | Right | 264 lbs | | 27.3 | 25.1 | MHFKR 195 |
| Mercedes | B Class | 2019 - 2025 | 5 | Right | 175 lbs | | 26.7 | N/A | MHFKR 436 |
| Mercedes | C Class | 2008 - 2014 | 5 | Right | 220 lbs | | | N/A | |
| Mercedes | C Class Wagon | 2015 - 2021 | 5 | Right | 220 lbs | | 23.7 | N/A | MHFKR 318 |
| Mercedes | E Class Estate | 2017 - 2023 | 5 | Right | 264 lbs | | 24 | 23 | MHFKR 353 |
| Mercedes | GL Class | 2007 - 2012 | 7 | Right | 264 lbs | | 30.7 | 28.7 | |
| Mercedes | GLA Class | 2020 - 2026 | 5 | Right | 220 lbs | | 26.3 | N/A | MHFKR 499 |
| Mercedes | GLB | 2020 - 2026 | 5 | Right | 220 lbs | 30.1 | 27.9 | N/A | MHFKR 584 |
| Mercedes | GLC Class | 2016 - 2022 | 5 | Left | 175 lbs | | 27.4 | N/A | MHFKR 424 |
| Mercedes | GLE Class | 2015 - 2018 | 5 | Right | 264 lbs | | 29.1 | 28.3 | MHFKR 356 |
| Mercedes | GLE Class | 2019 - 2026 | 5 | Right | 352 lbs | | 29.3 | 26.9 | MHFKR 437 |
| Mercedes | GLS | 2016 -2019 | 7 | Right | 440 lbs | | 30.5 | 28.3 | MHFKR 382 |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Mercedes | GLS | 2020 - 2026 | 5 | Left | 440 lbs | | 30.3 | 27.1 | |
| Mercedes | M Class | 2001 - 2005 | 5 | Right | 220 lbs | | | N/A | |
| Mercedes | M Class | 2006 - 2011 | 5 | Right | 330 lbs | | 31.1 | 29.1 | MHFKR 188 |
| Mercedes | M Class | 2012 - 2015 | 5 | Right | 352 lbs | | 30.3 | 28.3 | MHFKR 201 |
| Mercedes | Metris | 2016 - 2023 | 8 | Right | 440 lbs | | 32.6 | 30.6 | MHFKR 323 |
| MINI | Countryman | 2017 - Mid 2020 | 5 | Right | 175 lbs | | 25.7 | N/A | MHFKR 381 |
| Mitsubishi | Eclipse Cross | 2018 - 2022 | 5 | Right | 175 lbs | | 26.1 | N/A | MHFKR 388 |
| Mitsubishi | Outlander | 2014 - 2021 | 7 | Right | 330 lbs | | 28 | 25.9 | MHFKR 209 |
| Mitsubishi | Outlander PHEV | 2014 - 2021 | 5 | Right | 330 lbs | | 28 | 25.9 | |
| Mitsubishi | Outlander PHEV | 2022 - 2026 | 7 | Right | 175 lbs | N/A | 28.7 | N/A | |
| Mitsubishi | Outlander Sport | 2011 - 2023 | 5 | Right | 175 lbs | 24.8 | 24.6 | | |
| Nissan | Cube | 2009 - 2014 | 5 | Left | 220 lbs | | 31.9 | N/A | |
| Nissan | Kicks | 2018 - 2024 | 5 | Left | 175 lbs | 28.5 | 25.8 | N/A | |
| Nissan | Murano | 2003 - 2007 | 5 | Left | 440 lbs | | 28.5 | 26.7 | MHFKR 401 |
| Nissan | Murano | 2015 - 2024 | 5 | Left | 220 lbs | 28.3 | 26.7 | N/A | |
| Nissan | NV200 | 2010 - 2020 | 7 | Right | 330 lbs | | 31.5 | 29.5 | |
| Nissan | Pathfinder | 2005 - 2012 | 7 | Right | 220 lbs | | 31.5 | N/A | |
| Nissan | Pathfinder | 2013 - 2021 | 7 | Right | 440 lbs | | 29.9 | 27.9 | |
| Nissan | Pathfinder | 2022 - 2026 | 7 | Right | 330 lbs | | 29.1 | 26.3 | |
| Nissan | Quest | 2011 - 2017 | 7 | Right | 440 lbs | | 34 | 32.5 | |
| Nissan | Rogue | 2014 - 2020 | 5 | Right | 330 lbs | | 28.9 | 27.5 | MHFKR 274 |
| Nissan | Rogue | 2014 - 2020 | 7 | Right | 330 lbs | | 28.8 | 26.7 | MHFKR 279 |
| Nissan | Rouge | 2021 - 2026 | 5 | Right | 275 lbs | | 29.1 | 27.5 | |
| Nissan | Rouge Sport/Qashqai | 2017 - 2022 | 5 | Right | 220 lbs | | 23.6 | N/A | MHFKR 262 |
| Nissan | Xterra | 2005 - 2012 | 5 | Right | 440 lbs | | 31.9 | 29.9 | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Nissan | Versa Note | 2014 - 2019 | 5 | Right | 220 lbs | | 25.6 | N/A | MHFKR 240 |
| Porsche | Cayenne | 2011 - 2018 | 5 | Right | 275 lbs | | 27.6 | 25.6 | MHFKR 255 |
| Porsche | Cayenne | 2019 - 2026 | 5 | Right | 275 lbs | | 27.6 | 25.6 | MHFKR 255 |
| Ram | Promaster City | 2015 - 2025 | 5 & 7 | Right | 330 lbs | | 30 | 27.9 | |
| Subaru | Ascent | 2019 - 2026 | 7 | Right | 275 lbs | | 29.3 | 26.7 | |
| Subaru | Crosstrek | 2018 - 2023 | 5 | Left | 220 lbs | 25.6 | 23.2 | | |
| Subaru | Forester | 2003 - 2008 | 5 | Right | 220 lbs | | 28 | N/A | |
| Subaru | Forester | 2014 - 2018 | 5 | Left | 330 lbs | | 29.6 | 27.5 | MHFKR 333 |
| Subaru | Forester | 2019 - 2024 | 5 | Left | 275 lbs | | 29.7 | 27.9 | |
| Subaru | Forester | 2025 - 2028 | 5 | Left | 220 lbs | N/A | 29.5 | N/A | |
| Subaru | Forester XT | 2012 - 2018 | 5 | Left | 330 lbs | | 29.2 | 27.1 | |
| Subaru | Outback | 2005 - 2009 | 5 | Right | 220 lbs | | 26.4 | N/A | |
| Subaru | Outback | 2010 - 2014 | 5 | Left | 220 lbs | | 27.6 | N/A | MHFKR 327 |
| Subaru | Outback | 2015 - 2019 | 5 | Right | 330 lbs | | 27.5 | 25.5 | |
| Subaru | Outback | 2020 - 2025 | 5 | Left | 220 lbs | | 27.5 | N/A | |
| Tesla | Model X | 2016 - 2026 | 6 | Right | 175 lbs | | 23.8 | N/A | |
| Tesla | Model Y | 2019 - 2026 | 5 | Right | 175 lbs | 25 | n/a | N/A | |
| Toyota | 4 Runner | 2010 - 2024 | 5 | Left | 220 lbs | | 30.7 | | |
| Toyota | BZ4X | 2023 - 2026 | 5 | Left | 264 lbs | 25.2 | 25.7 | N/A | |
| Toyota | C-HR | 2018 - 2023 | 5 | Left | 175 lbs | 22.6 | n/a | N/A | MHFKR 379 |
| Toyota | Corolla Cross | 2022 - 2026 | 5 | Left | 175 lbs | 27.2 | 24 | N/A | |
| Toyota | Highlander | 2001 - 2007 | 5 | Right | 330 lbs | | 28.4 | 26.3 | |
| Toyota | Highlander | 2014 - 2019 | 7 | Left | 264 lbs | | 29.6 | 27.5 | |
| Toyota | Highlander | 2020 - 2026 | 7 | Left | 220 lbs | | 27.9 | N/A | |
| Toyota | Grand Highlander | 2024 - 2026 | 7 | Left | 275 lbs | 29.5 | 29.5 | 26.7 | |
| Toyota | Land Cruiser | 1999 - 2007 | 8 | Left | 220 lbs | | | N/A | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| Toyota | Land Cruiser | 2008 - 2021 | 7 | Left | 220 lbs | | 30.7 | N/A | MHFKR 200 |
| Toyota | Rav4 | 2013 - 2018 | 5 | Right | 352 lbs | | 28.4 | 26.3 | MHFKR 241 |
| Toyota | Rav4 (NA Version) | 2006 - 2012 | 5 | Left | 440 lbs | | 34.3 | 32.2 | MHFKR 211 |
| Toyota | Rav4 (NA Version) | 2019 - 2026 | 5 | Left | 440 lbs | | 31.1 | 28.3 | |
| Toyota | Sienna | 2004 - 2010 | 7 | Right | 352 lbs | | 35.1 | 33 | |
| Toyota | Sienna | 2011 - 2020 | 7 | Right | 352 lbs | | 34.3 | 32.2 | |
| Toyota | Sienna | 2021 - 2026 | 7 | Right | 440 lbs | | 34 | 31.5 | |
| Toyota | Venza | 2021 - 2024 | 5 | Left | 220 lbs | 27 | 24 | N/A | |
| Toyota | Yaris (Hatchback) | 2012 - 2018 | 5 | Right | 175 lbs | 25.6 | n/a | N/A | MHFKR 345 |
| Volvo | V50 | 2004 - 2011 | 5 | Right | 220 lbs | | 23.7 | N/A | |
| Volvo | V70 | 2001 - 2007 | 5 | Right | 220 lbs | | | N/A | |
| Volvo | V70 | 2008 - 2016 | 5 | Right | 264 lbs | | | | |
| Volvo | V90 | 2017 - 2021 | 5 | Right | 175 lbs | | 25.3 | N/A | MHFKR 443 |
| Volvo | XC40 | 2018 - 2026 | 5 | Left | 175 lbs | | 26.2 | N/A | MHFKR 414 |
| Volvo | XC60 | 2009 - 2017 | 5 | Right | 352 lbs | | 27.2 | 25.1 | |
| Volvo | XC60 | 2018 - 2026 | 5 | Left | 175 lbs | | 28.3 | N/A | MHFKR 386 |
| Volvo | XC70 | 2001 - 2007 | 5 | Right | 220 lbs | | | N/A | |
| Volvo | XC70 | 2008 - 2016 | 5 | Right | 264 lbs | | | | |
| Volvo | XC90 | 2015 - 2026 | 7 | Right | 440 lbs | | 27.5 | 25.9 | |
| VW | Atlas | 2018 - 2026 | 5 | Right | 275 lbs | | 30.3 | 27.7 | |
| VW | Atlas Cross Sport | 2020 - 2026 | 5 | Right | 220 lbs | | 26.1 | N/A | |
| VW | Golf Hatchback | 2015 - 2020 | 5 | Left | 175 lbs | 23 | n/a | N/A | MHFKR 354 |
| VW | Golf Mk8(eTSi) | 2022 - 2025 | 5 | Left | 175 lbs | 25.9 | n/a | N/A | MHFKR 559 |
| VW | Golf Wagon | 2015 - 2025 | 5 | Right | 220 lbs | | 24.4 | N/A | MHFKR 256 |
| VW | Golf Wagon | 2010 - 2014 | 5 | Right | 330 lbs | | 27.2 | 25.1 | |
| VW | ID.4 | 2021 - 2026 | 5 | Left | 175 lbs | 28.7 | 24.6 | | |

Smart Lifter Vehicle Specific Kits *Please note Universal A Frame kits may be available*

| Vehicle Make | Vehicle Model | Vehicle Model Year | # of Vehicle Seats | Vehicle Installation Side | Max Mobility Device Weight | Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift | Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift | Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift | Optional Mobility Device Restraint Kit <small>*universal option available*</small> |
|---------------------|----------------------|---------------------------|---------------------------|----------------------------------|-----------------------------------|--|---|---|--|
| VW | Taos | 2022 - 2026 | 5 | Left | 275 lbs | | 28 | 25.6 | |
| VW | Tiguan (LWB) | 2018 - 2022 | 7 | Right | 220 lbs | | 25.9 | N/A | MHFKR 422 |
| VW | Tiguan (SWB) | 2018 - 2025 | 5 | Left | 220 lbs | | 27.3 | N/A | MHFKR 359 |
| VW | Touareg | 2011 - 2018 | 5 | Right | 352 lbs | | 28.4 | 26.3 | MHFKR 284 |