

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Acura	MDX	2007 - 2013	7	Right	220 lbs		29.6	N/A	
Acura	MDX	2014 - 2020	7	Right	220 lbs		27.9	N/A	
Acura	MDX	2022 - 2026	7	Right	220 lbs	29.5	27.5	N/A	
Acura	RDX	2019 - 2026	5	Right	440 lbs		27.5	24.8	
Audi	Q3	2015 - 2018	5	Left	175 lbs	24.8	n/a	N/A	MHFKR 264
Audi	Q3	2019 - 2025	5	Right	220 lbs	24.3	24.3	N/A	MHFKR 427
Audi	Q4 E-Tron	2022 - 2026	5	Right	175 lbs	25.7	23	N/A	
Audi	Q5	2009 - 2017	5	Right	275 lbs		30.3	26.3	MHFKR 227
Audi	Q5	2018 - 2025	5	Left	275 lbs		26.3	24.2	MHFKR 374
Audi	Q7	2016 - 2025	7	Right	175 lbs		27.3	N/A	MHFKR 461
BMW	3 Series (Wagon)	2012 - 2018	5	Right	175 lbs	23.8	n/a	N/A	MHFKR 339
BMW	X1	2012 - 2015	5	Right	220 lbs		24	22	MHFKR 189
BMW	X1	2016 - 2022	5	Right	175 lbs		25	N/A	MHFKR 342
BMW	X2	2019 - 2025	5	Right	175 lbs	26.5	n/a	N/A	MHFK 439
BMW	X3	2004 - 2010	5	Right	440 lbs		30.5	28.5	MHFKR 208
BMW	X3	2011 - 2017	5	Right	440 lbs		29.2	27.1	
BMW	X3	2018 - 2024	5	Left	175 lbs	28.9	25.3	N/A	MHFKR 448
BMW	X5	2019 - 2026	5	Right	275 lbs		29.3	27.1	
BMW	X7	2019 - 2026	7	Right	220 lbs	N/A	29.2	N/A	
Buick	Enclave	2008 - 2017	7	Right	440 lbs		28	25.9	
Buick	Enclave	2018 - 2024	7	Right	440 lbs		29.1	27.3	
Buick	Enclave	2025 - 2028	6	Right	220 lbs	N/A	26.3	N/A	
Buick	Encore	2013 - 2022	5	Right	220 lbs		26.8	24.8	MHFKR 252
Buick	Encore GX	2020 - 2026	5	Left	175 lbs	25.9	n/a	N/A	
Buick	Envision	2021 - 2026	5	Left	175 lbs		23.8	N/A	
Buick	Envista	2024 - 2028	5	Left	175 lbs	N/A	22.4	N/A	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Buick	Rendezvous	2002 - 2007	5	Right	275 lbs		33.5	31.4	
Buick	Envision	2016 - 2020	5	Right	275 lbs		27.9	25.6	
Cadillac	Lyric	2023 - 2027	5	Right	175 lbs	28.8	25.3	N/A	
Cadillac	SRX	2010 - 2016	5	Right	275 lbs		25.2	23.2	
Cadillac	XT5	2017 - 2026	5	Right	220 lbs		27.1	N/A	
Cadillac	XT6	2020 - 2025	7	Left	275 lbs		27.7	25.1	
Chevrolet	Blazer	2019 - 2025	5	Right	220 lbs		27.5	N/A	
Chevrolet	Blazer EV	2024 - 2027	5	Right	175 lbs	N/A	24	N/A	
Chevrolet	Bolt EUV	2022 - 2023	5	Left	175 lbs	n/a	22.8	N/A	
Chevrolet	Captiva	2007 - 2018	7	Right	220 lbs		30.4	N/A	
Chevrolet	City Express	2015 - 2018	7	Right	275 lbs		31.5	29.5	
Chevrolet	Equinox	2005 - 2009	5	Right	220 lbs		27.8	25.7	
Chevrolet	Equinox	2010 - 2017	5	Right	440 lbs		29.6	27.5	
Chevrolet	Equinox	2018 - 2024	5	Right	220 lbs		27.5	25	
Chevrolet	Equinox	2025 - 2028	5	Left	220 lbs	N/A	26.3	N/A	
Chevrolet	Express	1996 - 2026	12	Right	275 lbs		33.5	31.4	
Chevrolet	HHR	2006 - 2011	5	Right	275 lbs		26.8	24.8	
Chevrolet	Orlando	2011 - 2018	7	Right	220 lbs		27.2	N/A	
Chevrolet	Suburban	2007 - 2014	7	Right	440 lbs		34	32	
Chevrolet	Suburban	2021 - 2026	7	Right	275 lbs		29.5	28.3	
Chevrolet	Tahoe	2007 - 2014	7	Right	440 lbs		31.9	29.9	
Chevrolet	Tahoe	2015 - 2020	7	Right	440 lbs		27	25	
Chevrolet	Tahoe	2021 - 2026	5	Right	220 lbs		29.9	N/A	
Chevrolet	Trailblazer	2002 - 2009	5	Right	220 lbs		28.8	26.7	
Chevrolet	Trailblazer	2021 - 2026	5	Left	220 lbs	26.7	25.2	N/A	
Chevrolet	Traverse	2009 - 2017	7	Right	440 lbs		28	25.9	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Chevrolet	Traverse	2018 - 2023	7	Right	440 lbs		31.4	29.1	
Chevrolet	Traverse	2024 - Split Gen	6-Jul	Right	440 or 220	N/A	31.4 or 26.3	29.1 or N/A	
Chevrolet	Traverse	2025 - 2028	6	Right	220 lbs	N/A	26.3	N/A	
Chevrolet	Trax	2015 - 2023	5	Right	220 lbs		26.8	24.8	
Chevrolet	Trax	2024 - 2028	5	Left	175 lbs	N/A	21.6	N/A	
Chevrolet	Uplander	2005 - 2009	7	Right	440 lbs		35.1	33	
Chrysler	Pacifica	2017 - 2026	7	Right	440 lbs		34.6	30.5	
Chrysler	Town & Country w/ out Stow n Go	2001 - 2007	7	Right	275 lbs		32.3	30.3	
Chrysler	Town & Country w/ Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	
Chrysler	Town & Country w/ Stow n Go	2008 - 2016	7	Right	440 lbs		35.6	33.8	
Dodge	Caravan & Grand Caravan w/ Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	
Dodge	Caravan w/ out Stow n Go	2001 - 2007	7	Right	440 lbs		30.3	28.3	
Dodge	Durango	2004 - 2009	7	Right	275 lbs		29.6	27.5	
Dodge	Durango	2011 - 2026	6	Right	275 lbs	29.5	29.5	26.9	
Dodge	Grand Caravan w/ Stow n Go	2008 - 2020	7	Right	440 lbs		35.6	33.8	
Dodge	Journey	2009 -2020	5	Left	220 lbs		26.7	N/A	
Fiat	500L	2014 - 2022	5	Right	275 lbs		28.4	26.3	
Fiat	500X	2016 - 2022	5	Right	175 lbs	21	n/a	n/a	MHFKR 423
Ford	Bronco Sport	2021 - 2026	5	Right	220 lbs		26.5	N/A	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Ford	Ecosport	2018 - 2022	5	Right	220 lbs		30.9	N/A	MHFKR 396
Ford	Edge	2007 - 2014	5	Right	440 lbs	31.5	28.3	26	
Ford	Edge	2015 - 2024	5	Right	220 lbs		27.9	25.9	
Ford	Escape	2001 - 2007	5	Right	275 lbs		30.7	28.7	
Ford	Escape	2008 - 2012	5	Right	275 lbs		30.7	28.3	
Ford	Escape	2013 - 2019	5	Right	275 lbs		29.6	26.9	MHFKR 202
Ford	Escape	2020 - 2026	5	Right	220 lbs		27.5	N/A	MHFKR 477
Ford	E-Series	1992 - 2014	12	Right	440 lbs		37.4	35.4	
Ford	Expedition	2007 - 2017	7	Right	275 lbs		30.3	28.3	
Ford	Expedition	2018 - 2024	7	Right	220 lbs	34	29.9	N/A	
Ford	Expedition EL/Max	2007 - 2017	7	Right	440 lbs		31.1	29.1	
Ford	Expedition MAX	2018 - 2024	8	Right	220 lbs	31.8	29.1	N/A	MHFKR 674
Ford	Explorer	2002 - 2005	5	Left	440 lbs		30.7	28.7	
Ford	Explorer	2011 - 2019	7	Right	275 lbs		27.9	25.9	
Ford	Explorer	2020 - 2026	6	Left	440 lbs		30.5	27.5	
Ford	Flex	2013 - 2019	7	Right	440 lbs		29	27	
Ford	Focus Wagon	2012 - 2018	5	Right	275 lbs		28.2	26.1	
Ford	Freestar	2004 - 2007	7	Right	275 lbs		33.9	31.8	
Ford	Mustang Mach-E	2021 - 2025	5	Right	175 lbs	27.1	24		
Ford	Transit Connect LWB	2014 - 2023	7	Right	440 lbs		33.5	31.49	MHFKR 324
Ford	Transit Connect SWB	2014 - 2023	5	Right	440 lbs		33.5	31.4	
Genesis	GV80	2021 - 2026	5	Right	220 lbs		26.7	N/A	
GMC	Acadia	2007 - 2017	7	Right	440 lbs		28	25.9	
GMC	Acadia	2017 - 2023	7	Left	275 lbs		28.3	25.5	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

Vehicle Make	Vehicle Model	Vehicle Model Year	# of Vehicle Seats	Vehicle Installation Side	Max Mobility Device Weight	Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift	Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift	Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift	Optional Mobility Device Restraint Kit <small>*universal option available*</small>
GMC	Acadia	2024 - 2028	6	Right	220 lbs	N/A	26.3	N/A	
GMC	Envoy	2002 - 2009	5	Right	220 lbs		28.8	26.7	
GMC	Savana	1996 - 2002	12	Right	275 lbs		33.5	31.4	
GMC	Terrain	2010 - 2017	5	Right	440 lbs		29.6	27.5	
GMC	Terrain	2018 - 2024	5	Right	220 lbs		27.5	25	
GMC	Terrain	2025 - 2027	5	Left	220 lbs	N/A	26.3	N/A	
GMC	Yukon	2015 - 2020	7	Right	440 lbs		27	25	
GMC	Yukon	2021 - 2026	5	Right	220 lbs		29.9	N/A	
GMC	Yukon Denali	2007 - 2014	7	Right	440 lbs		31.9	29.9	
GMC	Yukon XL	2007 - 2014	7	Right	440 lbs		34	32	
Honda	CR-V	2002 - 2006	5	Left	220 lbs			N/A	
Honda	CR-V	2007 - 2011	5	Right	220 lbs		31.5	N/A	
Honda	CR-V	2012 - 2016	5	Right	220 lbs		30.3	28.3	MHFKR 203
Honda	CR-V	2017 - 2022	5	Right	440 lbs		30.9	28.7	
Honda	CR-V	2023 - 2026	5	Right	440 lbs	31.1	29.8	27.6	
Honda	Fit	2015 - 2020	5	Right	220 lbs		27.6	N/A	MHFKR 321
Honda	HR-V	2016 - 2022	5	Right	440 lbs		26.6	24.6	MHFKR 313
Honda	HR-V	2023 - 2027	5	Right	220 lbs	30.7	26.9		
Honda	Odyssey	2005 - 2010	7	Right	275 lbs		35.9	33.8	
Honda	Odyssey	2011 - 2013	7	Right	275 lbs		35.5	33.4	
Honda	Odyssey	2014 - 2017	7	Right	440 lbs		37	35	
Honda	Odyssey	2018 - 2026	7	Right	440 lbs		33.4	31.4	
Honda	Passport	2019 - 2025	5	Right	440 lbs		29.1	27.1	
Honda	Pilot	2009 - 2015	7	Right	440 lbs		29.6	27.5	
Honda	Pilot	2016 - 2022	8	Right	275 lbs		29.1	26.3	
Honda	Pilot	2023 - 2026	7	Right	275 lbs	29.1	28.3	25.9	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

Vehicle Make	Vehicle Model	Vehicle Model Year	# of Vehicle Seats	Vehicle Installation Side	Max Mobility Device Weight	Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift	Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift	Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift	Optional Mobility Device Restraint Kit <small>*universal option available*</small>
Honda	Fit	2009 - 2014	5	Right	220 lbs			N/A	
Hyundai	Entourage	2007 - 2009	7	Left	220 lbs			N/A	
Hyundai	Entourage w/ Stow n Go	2006 - 2009	7	Right	275 lbs		30.7	28.7	
Hyundai	Ioniq 5	2022 - 2026	5	Left	175 lbs	25.8	N/A	N/A	
Hyundai	Kona (Gas/HEV/BEV)	2018 - 2023	5	Right	175 lbs	27	25	N/A	MHFKR 431
Hyundai	Kona (mHEV)	2019 - 2023	5	Right	175 lbs	27	25	N/A	
Hyundai	Kona (Gas)	2024 - 2026	5	Right	220 lbs	28.3	26.5	N/A	
Hyundai	Palisade	2020 - 2025	7	Right	440 lbs		29.5	26.7	
Hyundai	Palisade	2026 - 2027	7	Right	220 lbs	31.25	29.5	N/A	
Hyundai	Santa Fe	2001 - 2006	5	Right	220 lbs			N/A	
Hyundai	Santa Fe	2007 - 2012	5	Right	220 lbs		31.5	N/A	
Hyundai	Santa Fe	2019 - 2020	5	Right	275 lbs		27.5	25.1	
Hyundai	Santa Fe	2021 - 2023	5	Right	275 lbs		27.5	26.7	
Hyundai	Santa Fe	2024 - 2028	7	Left	440 lbs	30.7	29.3	27.3	
Hyundai	Santa Fe Sport	2013 - 2018	5	Right	220 lbs		27.8	25.7	
Hyundai	Tucson	2005 - 2009	5	Left	275 lbs		28.4	26.3	
Hyundai	Tucson	2010 - 2015	5	Right	275 lbs				
Hyundai	Tucson	2016 - 2021	5	Left	220 lbs		26.7	N/A	
Hyundai	Tucson HEV	2022 - 2026	5	Left	220 lbs	29.5	27.5	N/A	
Infiniti	QX60	2022 - 2026	8	Right	220 lbs	30.0	28.0	N/A	
Jaguar	E-Pace	2018 - 2026	5	Right	175 lbs		24.6	N/A	MHFKR 393
Jaguar	F Pace	2016 - 2026	5	Right	175 lbs		26.1	N/A	MHFKR 361
Jeep	Cherokee	2014 - 2023	5	Right	220 lbs		25.7	N/A	
Jeep	Compass	2018 - 2025	5	Right	175 lbs		24.6	N/A	MHFKR 395

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Jeep	Grand Cherokee	2005 - 2010	5	Right	220 lbs			N/A	
Jeep	Grand Cherokee	2011 - 2020	5	Right	220 lbs		28	25.9	MHFKR 254
Jeep	Grand Cherokee	2022 - 2026	5	Right	275 lbs		28.9	26.5	
Jeep	Grand Cherokee L	2021 - 2026	6	Right	275 lbs		30.3	27.2	
Jeep	Liberty	2008 - 2012	5	Right	220 lbs		26.6	24.6	
Jeep	Patriot	2007 - 2017	5	Left	175 lbs	n/a	23.2	N/A	
Jeep	Renegade	2015 - 2020	5	Right	220 lbs		27	25	
Jeep	Wagoneer/Grand Wagoneer	2022 - 2026	8	Left	275 lbs		31.5	30.5	
Jeep	Wrangler	2018 - 2026	5	Right	175 lbs	30.2	27.5	N/A	
Kia	Borrego	2009 - 2011	7	Right	440 lbs		30.9	28	
Kia	Carnival	2022 - 2026	7	Right	440 lbs		35.2	31.8	
Kia	Niro HEV	2017 - 2022	5	Left	175 lbs		22.4	N/A	
Kia	Rondo/Carens	2006 - 2012	5 or 7	Right	220 lbs				
Kia	Rondo/Carens	2013 - 2017	5	Right	220 lbs		27.6	N/A	MHFKR 233
Kia	Sedona	2002 - 2005	7	Left	220 lbs			N/A	
Kia	Sedona NON Stow n Go	2006 - 2014	7	Left	220 lbs			N/A	
Kia	Sedona with Stow n Go	2006 - 2014	7	Right	275 lbs		30.7	28.7	
Kia	Seltos	2021 - 2026	5	Right	220 lbs		25.9	N/A	
Kia	Sorento	2003 - 2009	5	Right	220 lbs		30	N/A	
Kia	Sorento	2011 - 2015	7	Right	275 lbs		28.4	26.3	
Kia	Sorento	2016 - 2020	7	Left	275 lbs		27.8	25.7	MHFKR 319
Kia	Sorento	2021 - 2026	7	Left	440 lbs		27.6	25.4	MHFKR 537
Kia	Soul	2010 - 2013	5	Right	220 lbs		26	N/A	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Kia	Soul	2014 - 2019	5	Right	220 lbs		25.8	N/A	MHFKR 331
Kia	Soul	2020 - 2026	5	Right	220 lbs		25.9	N/A	
Kia	Sportage	2011 - 2016	5	Right	175 lbs			N/A	MHFKR 192
Kia	Sportage	2017 - 2022	5	Left	175 lbs		25.9	N/A	MHFKR 332
Kia	Sportage LWB	2023 - 2026	5	Left	275 lbs		29.1	27.1	
Kia	Sportage HEV LWB	2023 - 2026	5	Left	275 lbs	N/A	28.3	25.9	
Kia	Telluride	2020 - 2026	8	Right	440 lbs		29.5	27.9	
Kia	Telluride	2027 - 2028	7	Right	220 lbs	29.5	27.5	N/A	
Kia	Sedona	2015 - 2021	7	Right	440 lbs		33	27.5	
Land Rover	Discovery 2	1999 - 2003	5	Left	220 lbs			N/A	
Land Rover	Discovery 5	2017 - 2026	7	Right	220 lbs		30.7	N/A	MHFKR 363
Land Rover	Discovery Sport	2014 - 2019	7	Right	440 lbs		28.6	26.5	MHFKR 347
Land Rover	Freelander	2000 - 2005	5	Left	220 lbs			N/A	
Land Rover	Freelander	2006 - 2014	5	Right	275 lbs		28.2	26.1	
Land Rover	Range Rover Evoque	2011 - 2018	5	Right	220 lbs		24.8	N/A	MHFKR 190
Land Rover	Range Rover Evoque	2019 - 2023	5	Right	175 lbs	29.1	24.2	N/A	
Land Rover	Range Rover Sport	2005 - 2013	5	Right	220 lbs		28.6	N/A	
Land Rover	Range Rover Sport	2014 - 2022	5	Right	220 lbs		28	25.9	MHFKR 249
Land Rover	Range Rover Sport	2023-2027	5	Right	220 lbs		28	25.9	MHFKR 645
Land Rover	Range Rover Velar	2018 - 2026	5	Right	175 lbs		23	N/A	MHFKR 487
Lexus	NX	2015 - 2021	5	Right	175 lbs		24.2	N/A	MHFKR 419
Lexus	NX	2022 - 2026	5	Right	220 lbs	27.6	24.5	N/A	
Lexus	RX	2010 - 2015	5	Right	220 lbs		25.4	24.4	MHFKR 250
Lexus	RX	2016 - 2022	5	Left	175 lbs		24	N/A	MHFKR 334

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Lexus	RX	2023 - 2026	5	Left	220 lbs	29.1	26.3	N/A	
Lexus	RX-L	2017 - 2022	5	Right	175 lbs	28	25.2	N/A	
Lincoln	Aviator	2020 - 2026	7	Left	220 lbs		27.1	N/A	
Lincoln	Corsair	2020 - 2026	5	Right	220 lbs		25.5	N/A	
Lincoln	MKT	2010 - 2019	6	Right	220 lbs	N/A	27.5	N/A	
Lincoln	MKX	2016 - 2018	5	Right	275 lbs		25.5	24.4	
Lincoln	Nautilus	2019 - 2023	5	Right	220 lbs	28.7	25.5	N/A	
Lincoln	Nautilus	2024 - 2026	5	Right	220 lbs	N/A	25.9	N/A	
Lincoln	Navigator	2025 - 2028	7	Right	220 lbs	N/A	29.9	N/A	
Mazda	5	2012 - 2015	7	Right	440 lbs		29.4	27.3	
Mazda	6	2003 - 2008	5	Right	220 lbs			N/A	
Mazda	CX-30	2020 - 2026	5	Right	175 lbs		22.6	N/A	MHFKR 475
Mazda	CX-50	2023 - 2026	5	Right	175 lbs	28.7	26.3	N/A	
Mazda	CX-70	2025 - 2028	5	Right	220 lbs	29.5	26.7	N/A	
Mazda	CX-5	2013 - 2016	5	Right	275 lbs		30.3	28.3	MHFKR 228
Mazda	CX-5	2017 - 2025	5	Right	275 lbs		27.5	25.1	
Mazda	CX-9	2007 - 2015	7	Right	275 lbs		29.6	27.5	
Mazda	MX-30 EV	2022 - 2026	5	Right	175 lbs	25.6	n/a	N/A	MHFKR 564
Mercedes	B Class	2006 - 2011	5	Right	220 lbs		27.6	N/A	
Mercedes	B Class	2012 - 2018	5	Right	220 lbs		27.3	25.1	MHFKR 195
Mercedes	B Class	2019 - 2025	5	Right	175 lbs		26.7	N/A	MHFKR 436
Mercedes	C Class	2008 - 2014	5	Right	220 lbs			N/A	
Mercedes	C Class Wagon	2015 - 2021	5	Right	220 lbs		23.7	N/A	MHFKR 318
Mercedes	E Class Estate	2017 - 2023	5	Right	220 lbs		24	23	MHFKR 353
Mercedes	GL Class	2007 - 2012	7	Right	220 lbs		30.7	28.7	
Mercedes	GLA Class	2020 - 2026	5	Right	220 lbs		26.3	N/A	MHFKR 499

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Mercedes	GLB	2020 - 2026	5	Right	220 lbs	30.1	27.9	N/A	MHFKR 584
Mercedes	GLC Class	2016 - 2022	5	Left	175 lbs		27.4	N/A	MHFKR 424
Mercedes	GLE Class	2015 - 2018	5	Right	220 lbs		29.1	28.3	MHFKR 356
Mercedes	GLE Class	2019 - 2026	5	Right	275 lbs		29.3	26.9	MHFKR 437
Mercedes	GLS	2016 -2019	7	Right	440 lbs		30.5	28.3	MHFKR 382
Mercedes	GLS	2020 - 2026	5	Left	440 lbs		30.3	27.1	
Mercedes	M Class	2001 - 2005	5	Right	220 lbs			N/A	
Mercedes	M Class	2006 - 2011	5	Right	275 lbs		31.1	29.1	MHFKR 188
Mercedes	M Class	2012 - 2015	5	Right	275 lbs		30.3	28.3	MHFKR 201
Mercedes	Metris	2016 - 2023	8	Right	440 lbs		32.6	30.6	MHFKR 323
MINI	Countryman	2017 - Mid 2020	5	Right	175 lbs		25.7	N/A	MHFKR 381
Mitsubishi	Eclipse Cross	2018 - 2022	5	Right	175 lbs		26.1	N/A	MHFKR 388
Mitsubishi	Outlander	2014 - 2021	7	Right	275 lbs		28	25.9	MHFKR 209
Mitsubishi	Outlander PHEV	2014 - 2021	5	Right	275 lbs		28	25.9	
Mitsubishi	Outlander PHEV	2022 - 2026	5	Right	220 lbs	N/A	28.7	N/A	
Mitsubishi	Outlander Sport	2011 - 2023	5	Right	175 lbs	24.8	24.6		
Nissan	Cube	2009 - 2014	5	Left	220 lbs		31.9	N/A	
Nissan	Kicks	2018 - 2024	5	Left	175 lbs	28.5	25.8	N/A	
Nissan	Murano	2003 - 2007	5	Left	440 lbs		28.5	26.7	MHFKR 401
Nissan	Murano	2015 - 2024	5	Left	220 lbs	28.3	26.7	N/A	
Nissan	NV200	2010 - 2020	7	Right	275 lbs		31.5	29.5	
Nissan	Pathfinder	2005 - 2012	7	Right	220 lbs		31.5	N/A	
Nissan	Pathfinder	2013 - 2021	7	Right	440 lbs		29.9	27.9	
Nissan	Pathfinder	2022 - 2026	7	Right	275 lbs		29.1	26.3	
Nissan	Quest	2011 - 2017	7	Right	440 lbs		34	32.5	
Nissan	Rogue	2014 - 2020	5	Right	275 lbs		28.9	27.5	MHFKR 274

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Nissan	Rogue	2014 - 2020	7	Right	275 lbs		28.8	26.7	MHFKR 279
Nissan	Rouge	2021 - 2026	5	Right	275 lbs		29.1	27.5	
Nissan	Rouge Sport/Qashqai	2017 - 2022	5	Right	220 lbs		23.6	N/A	MHFKR 262
Nissan	Xterra	2005 - 2012	5	Right	440 lbs		31.9	29.9	
Nissan	Versa Note	2014 - 2019	5	Right	220 lbs		25.6	N/A	MHFKR 240
Porsche	Cayenne	2011 - 2018	5	Right	275 lbs		27.6	25.6	MHFKR 255
Porsche	Cayenne	2019 - 2026	5	Right	275 lbs		27.6	25.6	MHFKR 255
Ram	Promaster City	2015 - 2025	5 & 7	Right	275 lbs		30	27.9	
Subaru	Ascent	2019 - 2026	7	Right	275 lbs		29.3	26.7	
Subaru	Crosstrek	2018 - 2023	5	Left	220 lbs	25.6	23.2		
Subaru	Forester	2003 - 2008	5	Right	220 lbs		28	N/A	
Subaru	Forester	2014 - 2018	5	Left	275 lbs		29.6	27.5	MHFKR 333
Subaru	Forester	2019 - 2024	5	Left	275 lbs		29.7	27.9	
Subaru	Forester	2025 - 2028	5	Left	220 lbs	N/A	29.5	N/A	
Subaru	Forester XT	2012 - 2018	5	Left	275 lbs		29.2	27.1	
Subaru	Outback	2005 - 2009	5	Right	220 lbs		26.4	N/A	
Subaru	Outback	2010 - 2014	5	Left	220 lbs		27.6	N/A	MHFKR 327
Subaru	Outback	2015 - 2019	5	Right	275 lbs		27.5	25.5	
Subaru	Outback	2020 - 2025	5	Left	220 lbs		27.5	N/A	
Tesla	Model X	2016 - 2025	6	Right	175 lbs		23.8	N/A	
Tesla	Model Y	19 - 25 1st gen	5	Right	175 lbs	25	n/a	N/A	
Toyota	4 Runner	2010 - 2024	5	Left	220 lbs		30.7		
Toyota	BZ4X	2023 - 2026	5	Left	220 lbs	25.2	25.7	N/A	
Toyota	C-HR	2018 - 2023	5	Left	175 lbs	22.6	n/a	N/A	MHFKR 379
Toyota	Corolla Cross	2022 - 2026	5	Left	175 lbs	27.2	24	N/A	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
Toyota	Highlander	2001 - 2007	5	Right	275 lbs		28.4	26.3	
Toyota	Highlander	2014 - 2019	7	Left	220 lbs		29.6	27.5	
Toyota	Highlander	2020 - 2026	7	Left	220 lbs		27.9	N/A	
Toyota	Grand Highlander	2024 - 2026	7	Left	275 lbs	29.5	29.5	26.7	
Toyota	Land Cruiser	1999 - 2007	8	Left	220 lbs			N/A	
Toyota	Land Cruiser	2008 - 2021	7	Left	220 lbs		30.7	N/A	MHFKR 200
Toyota	Prius	2016 - 2022	5	Left	175 lbs	N/A	22.8	N/A	
Toyota	Rav4	2013 - 2018	5	Right	275 lbs		28.4	26.3	MHFKR 241
Toyota	Rav4 (NA Version)	2006 - 2012	5	Left	440 lbs		34.3	32.2	MHFKR 211
Toyota	Rav4 (NA Version)	2019 - 2025	5	Left	440 lbs		31.1	28.3	
Toyota	Rav4 (NA Version)	2026 - 2028	5	Left	220 lbs	32.2	29.5	N/A	
Toyota	Sienna	2004 - 2010	7	Right	275 lbs		35.1	33	
Toyota	Sienna	2011 - 2020	7	Right	275 lbs		34.3	32.2	
Toyota	Sienna	2021 - 2026	7	Right	440 lbs		34	31.5	
Toyota	Venza	2021 - 2024	5	Left	220 lbs	27	24	N/A	
Toyota	Yaris (Hatchback)	2012 - 2018	5	Right	175 lbs	25.6	n/a	N/A	MHFKR 345
Volvo	V50	2004 - 2011	5	Right	220 lbs		23.7	N/A	
Volvo	V70	2001 - 2007	5	Right	220 lbs			N/A	
Volvo	V70	2008 - 2016	5	Right	220 lbs				
Volvo	V90	2017 - 2021	5	Right	175 lbs		25.3	N/A	MHFKR 443
Volvo	XC40	2018 - 2026	5	Left	175 lbs		26.2	N/A	MHFKR 414
Volvo	XC60	2009 - 2017	5	Right	275 lbs		27.2	25.1	
Volvo	XC60	2018 - 2026	5	Left	175 lbs		28.3	N/A	MHFKR 386
Volvo	XC70	2001 - 2007	5	Right	220 lbs			N/A	
Volvo	XC70	2008 - 2016	5	Right	220 lbs				
Volvo	XC90	2015 - 2026	7	Right	440 lbs		27.5	25.9	

**Smart Lifter Vehicle Specific Kits** \*Please note Universal A Frame kits may be available\*

<b>Vehicle Make</b>	<b>Vehicle Model</b>	<b>Vehicle Model Year</b>	<b># of Vehicle Seats</b>	<b>Vehicle Installation Side</b>	<b>Max Mobility Device Weight</b>	<b>Max Mobility Device Height in inches - LM 80 lb or LM 175 lb Lift</b>	<b>Max Mobility Device Height in inches - LC 175 lb or LC 220 lb Lift</b>	<b>Max Mobility Device Height in inches - LP 275 lb or LP 440 lb Lift</b>	<b>Optional Mobility Device Restraint Kit</b> <small>*universal option available*</small>
VW	Atlas	2018 - 2026	5	Right	275 lbs		30.3	27.7	
VW	Atlas Cross Sport	2020 - 2026	5	Right	220 lbs		26.1	N/A	
VW	Golf Hatchback	2015 - 2020	5	Left	175 lbs	23	n/a	N/A	MHFKR 354
VW	Golf Mk8(eTSi)	2022 - 2025	5	Left	175 lbs	25.9	n/a	N/A	MHFKR 559
VW	Golf Wagon	2015 - 2025	5	Right	220 lbs		24.4	N/A	MHFKR 256
VW	Golf Wagon	2010 - 2014	5	Right	275 lbs		27.2	25.1	
VW	ID.4	2021 - 2026	5	Left	175 lbs	28.7	24.6		
VW	Taos	2022 - 2026	5	Left	275 lbs		28	25.6	
VW	Tiguan (LWB)	2018 - 2022	7	Right	220 lbs		25.9	N/A	MHFKR 422
VW	Tiguan (LWB)	2018 - 2024	5	Right	220 lbs	30.3	27.5	N/A	MHFKR 694
VW	Touareg	2011 - 2018	5	Right	275 lbs		28.4	26.3	MHFKR 284